

Transformation and Storage Solution (TSS)

Nicolas Schmelling

Abstract

Recipe for TSS buffer. Buffer can be used for TSS transformation of non-competent *E. coli* cells.

Citation: Nicolas Schmelling Transformation and Storage Solution (TSS). **protocols.io**

dx.doi.org/10.17504/protocols.io.gt3bwqn

Published: 19 Dec 2016

Materials

- ✓ Polyethylene Glycol 3350 by Contributed by users
- Luria-Bertani (LB) broth, makes 1L K488 by [Amresco](#)
- dimethylsulfoxide (DMSO) by [Sigma Aldrich](#)

Protocol

Buffer ingredients

Step 1.

Chemical	Volume/Mass	Final Conc.
PEG 8000 or PEG 3350	20 g	20% (w/v)
DMSO	10 ml	10% (w/v)
MgCl ₂ * H ₂ O	2.03 g	100 mM
LB Medium	add 100 ml (pH 6.5)	



REAGENTS

- ✓ Polyethylene Glycol 3350 by Contributed by users

Buffer preparation

Step 2.

Prepare 2x TSS without DMSO and autoclave.

Buffer preparation

Step 3.

Add sterile DMSO after autoclaving.



REAGENTS

dimethylsulfoxide (DMSO) by [Sigma Aldrich](#)

Buffer storage

Step 4.

Divide in 10 ml aliquots and store at 4 °C