

# M0371 no RE mixture

## New England Biolabs

### Abstract

**Citation:** New England Biolabs M0371 no RE mixture. **protocols.io**

dx.doi.org/10.17504/protocols.io.crmv45

**Published:** 03 Feb 2015

## Protocol

### Step 1.

DNA - **1 pmol** of DNA ends

### Step 2.

CutSmart Buffer (10X), **2 µl**

 **AMOUNT**

2 µl Additional info:

### Step 3.

rSAP, **1 unit**

### Step 4.

purified H<sub>2</sub>O **to 20 µl**