

Iodoacetamide Solution

Kristin Corrier/Nathan VerBerkmoes

Abstract

For use in "[FASP Kit Protocol-ORNL Developed for Bacteriophage](#)"

Citation: Kristin Corrier/Nathan VerBerkmoes Iodoacetamide Solution. **protocols.io**
dx.doi.org/10.17504/protocols.io.dej3cm

Published: 21 Jan 2016

Protocol

Step 1.

100µL Urea solution (no DTT, no SDS) added to 1 tube iodoacetamide (provided with FASP kit). Pipette up and down 10-15 times to mix well and dissolve.