

TBST and 5% Milk

Alan Cone

Abstract

Citation: Alan Cone TBST and 5% Milk. [protocols.io](https://doi.org/10.17504/protocols.io.fw3bpqn)

[dx.doi.org/10.17504/protocols.io.fw3bpqn](https://doi.org/10.17504/protocols.io.fw3bpqn)

Published: 27 Sep 2016

Protocol

Step 1.

20 mM Tris pH 7.5

Step 2.

150 mM NaCl

Step 3.

0.1% Tween 2

Step 4.

5% Dry Milk