

Destain Buffer

Alan Cone

Abstract

Citation: Alan Cone Destain Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.fw9bph6

Published: 27 Sep 2016

Protocol

Step 1.

To a 500 mL Bottle, add 350 mL Water

 **AMOUNT**

350 ml Additional info:

 **REAGENTS**

✓ Distilled Water by Contributed by users

Step 2.

Add 100 mL 100% Methanol

 **AMOUNT**

100 ml Additional info:

 **REAGENTS**

 Methanol [PA-33900HPLCCS4L](#) by [P212121](#)

Step 3.

Add 50 mL Glacial Acetic Acid

 **AMOUNT**

50 ml Additional info:

 **REAGENTS**

Acetic acid, glacial [537020](#) by [Sigma Aldrich](#)