

Urea Solution

Kristin Corrier/Nathan VerBerkmoes

Abstract

For use in '[FASP Kit Protocol-ORNL Developed for Bacteriophage](#)'

Citation: Kristin Corrier/Nathan VerBerkmoes Urea Solution. **protocols.io**
[dx.doi.org/10.17504/protocols.io.dei3cd](https://doi.org/10.17504/protocols.io.dei3cd)

Published: 21 Jan 2016

Protocol

Step 1.

1mL Tris CaCl₂ or Tris HCl added to one tube of urea. Vortex until all powder dissolves.