



Feb 27, 2019

Working

Differentiation of NPC into cortical neurons

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x9bfr2n

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis

IPSC CORTICAL
DIFFERENTIATION
022017.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [IPSC CORTICAL DIFFERENTIATION](#) collection.

This method should be performed using sterile technique.

MATERIALS TEXT

Please refer to the attached full manuscript for required materials.

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

BEFORE STARTING

IMPORTANT: To generate cortical neurons, NPCs must be at passage 4 or lower at the time of plating for terminal differentiation. NPCs that are beyond passage 4 may exhibit higher densities of astrocyte contamination or inefficient neuronal differentiation.

- 1 Ensure cells are plated at an appropriate density (e.g.: 150K/well in 12-well plate; 75- 50K/well for 48 well plates; 30K/well for 8-well chamber slides) in NIM on pre-coated PLO/laminin plates.
- 2 After 24:00:00 , replace with cortical neuron differentiation medium (Neurobasal medium, 1x B27, 20 ng/mL BDNF, 20 ng/mL GDNF, 0.5mM cAMP, 1% Glutamax, 1% penicillin/streptomycin)
- 3 Feed cells every 2-3 days for 30 days.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited