



Oct 15, 2019

## Sanger sequencing

[iGEM Dusseldorf<sup>1</sup>](#)<sup>1</sup>Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.79ihr4e](https://doi.org/10.17504/protocols.io.79ihr4e)iGEM Dusseldorf 

### MATERIALS TEXT

barcodes

primers

template

water

- 1 Label 1.5 mL tubes with barcodes from GATC Biotech
- 2 Add template (400 to 500 ng), 2,5 µl of 10µM primer and water (10 uL total volume)
- 3 Place tubes inside a small plastic bag and drop off in GATC Collection Box



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited