

NEB 3X SDS Sample Loading Buffer

New England Biolabs

Abstract

Citation: New England Biolabs NEB 3X SDS Sample Loading Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.crcv2v

Published: 02 Feb 2015

Protocol

Step 1.

187.5 mM Tris-HCl (pH 6.8)

Step 2.

6% (w/v) SDS

Step 3.

30% glycerol

Step 4.

150 mM DTT

Step 5.

0.03% (w/v) bromophenol blue

Step 6.

2% β -mercaptoethanol