

Vitamin Solution for ESAW Media for Marine Phytoplankton

Ashley Humphrey

Abstract

Stock *vitamin solution* for ESAW Media for Marine Phytoplankton

Citation: Ashley Humphrey Vitamin Solution for ESAW Media for Marine Phytoplankton. **protocols.io**

dx.doi.org/10.17504/protocols.io.gdcbs2w

Published: 15 Nov 2016

Protocol

Vitamin solution, amounts per one mL

Step 1.

Thiamine

 **AMOUNT**

100 mg Additional info:

Vitamin solution, amounts per one mL

Step 2.

Vitamin B₁₂

 **AMOUNT**

2 mg Additional info:

Vitamin solution, amounts per one mL

Step 3.

Biotin

 **AMOUNT**

1 mg Additional info:

Vitamin solution

Step 4.

After mixing, filter sterilize with 0.2 µm filter and store in dark at -20°C. Divide into small aliquots if desired.