



Integrated analysis of DNA methylation profiling and gene expression profiling identifies novel markers in lung cancer in Xuanwei, China [↗](#)

PLOS One

Juan Wang¹

¹Department of Clinical Laboratory, the First Affiliated Hospital of Kunming Medical University, Kunming, China

[dx.doi.org/10.17504/protocols.io.rsxd6fn](https://doi.org/10.17504/protocols.io.rsxd6fn)



Juan Wang

Working
Oct 05, 2018

ABSTRACT

Background: Aberrant DNA methylation occurs frequently in cancer. The aim of this study was to identify novel methylation markers in lung cancer in Xuanwei, China, through integrated genome-wide DNA methylation and gene expression studies.

Methods: Differentially methylated regions (DMRs) and differentially expressed genes (DEGs) were detected on 10 paired lung cancer tissues and noncancerous lung tissues by methylated DNA immunoprecipitation combined with microarray (MeDIP-chip) and gene expression microarray analyses, respectively. Integrated analysis of DMRs and DEGs was performed to screen out candidate methylation-related genes. Both methylation and expression changes of the candidate genes were further validated and analyzed.

Results: Compared with normal lung tissues, lung cancer tissues expressed a total of 6,899 DMRs, including 5,788 hypermethylated regions and 1,111 hypomethylated regions. Integrated analysis of DMRs and DEGs identified 45 tumor-specific candidate genes: 38 genes whose DMRs were hypermethylated and expression was downregulated, and 7 genes whose DMRs were hypomethylated and expression was upregulated. The methylation and expression validation results identified 4 candidate genes (*STXBP6*, *BCL6B*, *FZD10*, and *HSPB6*) that were significantly hypermethylated and downregulated in most of the tumor tissues compared with the noncancerous lung tissues.

Conclusions: This integrated analysis of genome-wide DNA methylation and gene expression in lung cancer in Xuanwei revealed several genes regulated by promoter methylation that have not been described in lung cancer before. These results provide new insight into the carcinogenesis of lung cancer in Xuanwei and represent promising new diagnostic and therapeutic targets.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0203155>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Wang J, Duan Y, Meng Q, Gong R, Guo C, Zhao Y, Zhang Y (2018) Integrated analysis of DNA methylation profiling and gene expression profiling identifies novel markers in lung cancer in Xuanwei, China. PLoS ONE 13(10): e0203155. doi: [10.1371/journal.pone.0203155](https://doi.org/10.1371/journal.pone.0203155)

PROTOCOL STATUS

Working

GUIDELINES

The first part of the process requires quality control of experimental (IP) and control (input) samples to verify quality prior to microarray hybridization. You then independently label these samples using a NimbleGen Dual-Color DNA Labeling Kit and co-hybridize them to a NimbleGen DNA Methylation array using a NimbleGen Hybridization System. Following hybridization, the arrays are washed, dried, and scanned using a NimbleGen MS 200 Microarray Scanner. Array data are extracted and analyzed using our NimbleScan software and SignalMap software.

Then various verification experiments were carried out.

MATERIALS

NAME

CATALOG #













VENDOR



NimbleGen Arrays

Contributed by users

NAME 	CATALOG # 	VENDOR 
 NimbleGen HX3 Mixer (for 3x720K arrays; includes mixer port seals)	05 223 750 001	Contributed by users
 NimbleGen Dual-Color DNA Labeling Kit	06 370 250 001	Contributed by users
 NimbleGen Hybridization Kit	05 583 683 001	Contributed by users
 NimbleGen Sample Tracking Control Kit	05 223 512 001	Contributed by users
 NimbleGen Wash Buffer Kit	05 584 507 001	Contributed by users
 NimbleGen Array Processing Accessories	05 223 539 001	Contributed by users
 β-Mercaptoethanol	M3148	Contributed by users
 Absolute Ethanol 500 ml	E702-3	Contributed by users
 Isopropanol 500 ml	I-9516	Contributed by users
 Water: reagent grade, ACS, nonsterile, type 1	RC915025	Contributed by users
 QIAamp DNA Mini kit		Contributed by users
 QIAamp DNA Mini kit		Qiagen
 RNeasy Mini kit		Qiagen
 lysis of radio immunoprecipitation assay (RIPA)		Contributed by users
 phenylmethanesulfonyl fluoride (PMSF)		Contributed by users
 Transcriptor First Strand cDNA Synthesis kit		Contributed by users
 EpiTect Fast DNA Bisulfite kit		Contributed by users
 unmethylated: EpiTect Control DNA		Contributed by users
 methylated: EpiTect Control DNA		Qiagen
 EpiTect HRM PCR kit		Contributed by users
 Quantinova SYBR Green PCR Master Mix		Contributed by users
 sodium dodecyl sulfate–polyacrylamide gel		Contributed by users
 skimmed milk		Contributed by users
 TBST		Contributed by users
 Monoclonal rabbit anti-human STXBP6		Contributed by users
 Monoclonal rabbit anti-human BCL6B		Contributed by users
 Monoclonal rabbit anti-human FZD10		Contributed by users
 Monoclonal Mouse anti-human HSPB6		Contributed by users
 Monoclonal chicken anti-human GAPDH		Contributed by users
 Peroxidase-conjugated goat anti-rabbit IgG(HL)		Contributed by users
 Peroxidase-conjugated goat anti-mouse IgG(HL)		Contributed by users
 Peroxidase-conjugated goat anti-chicken IgG (HL)		Contributed by users

NAME	CATALOG #	VENDOR
 MassCLEAVE Kit	Part Number 10129	Contributed by users
 EZ DNA Methylation Kit		Contributed by users
 dd H ₂ O		Contributed by users
 10X Hot Star buffer		Qiagen
 dNTP mix		Contributed by users
 Hot Star taq		Qiagen
 Shrimp Alkaline Phosphatase		Contributed by users
 Polymerase buffer		Contributed by users
 T Cleavage Mix		Contributed by users
 DTT		Contributed by users
 T7 RNA & DNA Polymerase		Contributed by users
 RNase A		Contributed by users
 Clean Resin		Contributed by users

SAFETY WARNINGS

Wear gloves and take precautions to avoid sample contamination.

BEFORE STARTING

Cy dyes are light sensitive. Be sure to minimize light exposure of the dyes during use and store in the dark when not in use.

Cy dyes are ozone sensitive. Take the necessary precautions to keep atmospheric ozone levels below 20 ppb (parts per billion).

Cy dyes are humidity sensitive. Take the necessary precautions to keep humidity levels below 40%.

Roche NimbleGen has found that using VWR water and DTT for all post-hybridization washes results in higher signal from Cy dyes.

Roche NimbleGen recommends using a NanoDrop Spectrophotometer for quantifying and characterizing nucleic acid samples because this instrument requires only 1.5 µl of sample for analysis.

Perform all centrifugations at +15 to +25°C unless indicated otherwise.

MeDIP-chip

- Genomic DNA (gDNA) were extracted and purified from lung cancer and paired noncancerous lung tissues using the QIAamp DNA Mini kit (Qiagen, Hilden, Germany) according to the manufacturer's instructions.
For genome-wide methylation detection, equal quantities of gDNA from each of the paired tumor and normal lung specimens from 10 patients were pooled. A total of 6 µg of the pooled gDNA from each tissue type was digested by *MseI*, purified by using QIAquick PCR purification kit (Qiagen), and analyzed by agarose gel electrophoresis. *MseI*-digested gDNA (1.25 µg) was subjected to immunoprecipitation with a mouse monoclonal anti-5-methylcytosine antibody (Diagenode, Liège, Belgium) and then was purified through Qiagen MinElute columns and amplified by using a GenomePlex Complete Whole Genome Amplification (WGA2) kit from Sigma-Aldrich (St Louis, MO, USA). The amplified DNA samples were purified with a QIAquick PCR purification kit. The NimbleGen Dual-Color DNA Labeling kit (Roche NimbleGen Systems, Inc., Madison, WI, USA) was used for labeling according to the manufacturer's guidelines. Briefly, 1 µg DNA of each sample was incubated for 10 min at 98°C with 1 OD of Cy5-9mer primer (immunoprecipitation sample) or Cy3-9mer primer (input sample). Then, 100 pmol of deoxynucleoside triphosphates and 100 U of the Klenow fragment (New England Biolabs, Ipswich, MA, USA) were added, and the mix was incubated at 37°C for 3 h. The labeled DNA was purified by isopropanol/ethanol precipitation and hybridized to NimbleGen Human Meth 720K CpG RfSq Prom according to the manufacturer's instructions.

You can see the attachment of NG_DNA_Methylation_Guide_v7p2 for the specific experimental process.

 [ng_dna_methylation_guide_v7p2.pdf](#)

 [hb-0329-004-1102728-hb-qiaamp-dna-mini-blood-mini-0516-ww.pdf](#)

Gene expression microarray

- Gene expression profiling was performed on the same 10 paired samples used in MeDIP-chip analysis by using the Agilent Oligo Microarray Kit 8×60 K (Agilent Technologies, Santa Clara, CA, USA) as described previously. Significantly differentially expressed genes (DEGs) were identified by using the mixed model analysis of variance with a false discovery rate (Benjamini–Hochberg test)–adjusted q -value of ≤ 0.05 and absolute fold-change values ≥ 2 . Hierarchical clustering was carried out using cluster 3.0. GO analysis and Pathway analysis were performed to determine the roles of DEGs.

Methylation-sensitive high-resolution melting assay (MS-HRM)

- Genomic DNA (gDNA) were extracted and purified from lung cancer and paired noncancerous lung tissues using the QIAamp DNA Mini kit (Qiagen, Hilden, Germany) according to the manufacturer's instructions. gDNA was bisulfite-modified and purified by using the EpiTect Fast DNA Bisulfite kit (Qiagen) according to the manufacturer's instructions, and its concentration was adjusted to 10 ng/L. Methylation-sensitive high-resolution melting assay (MS-HRM) was performed on the 45 paired lung cancer and normal lung tissues to validate the methylation status of the candidate genes identified by the integrated analysis. A range of standards was included to control for bias in the sensitivity of the detection: 0% (unmethylated: EpiTect Control DNA; Qiagen), 100% (methylated: EpiTect Control DNA), and 50% (equal mixture of both templates). MS-HRM primers were designed to amplify a short amplicon size of less than 200 bp (<http://www.urogene.org/cgi-bin/methprimer>) and synthesized by Sangon Biotech (Shanghai, China) (S2 Table). MS-HRM assay was performed according to the instructions for the EpiTect HRM PCR kit (Qiagen) on the RotorGene Q (Qiagen) in triplicate. A sample amplification curve between the 50% and 100% standard curves was defined as hypermethylation, and a sample amplification curve between the 0% and 50% standard curves was defined as hypomethylation. You can see the attachments for the specific experimental process.

 [en-epitect-fast-bisulfite-conversion-handbook.pdf](#)

 [en-epitect-hrm-pcr-handbook.pdf](#)

 [hb-0329-004-1102728-hb-qiaamp-dna-mini-blood-mini-0516-ww.pdf](#)

Sequenom MassARRAY platform

- Quantitative methylation analysis of DMRs in the candidate genes was performed on the 45 paired tumor and normal lung tissues by using the Sequenom MassARRAY platform (CapitalBio) as described previously. This platform employs RNA base-specific cleavage and matrix-assisted laser desorption/ionization time-of-flight mass spectrometry. PCR primers for MassARRAY were designed by using EpiDesigner (<http://www.epidesigner.com>). The spectra methylation ratios were obtained via EpiTyper software version 1.0 (Sequenom, San Diego, CA, USA). You can see the attachment for the specific experimental process.

 [sequenom massarray platform.pdf](#)

Real-time fluorescent quantitative PCR

- Total RNA were extracted and purified from lung cancer and paired noncancerous lung tissues using RNeasy Mini kit (Qiagen, Hilden, Germany) according to the manufacturer's instructions. cDNA was synthesized by RNA reverse-transcription using the Transcriptor First Strand cDNA Synthesis kit (Roche, Mannheim, Germany). Real-time fluorescent quantitative PCR (RT-qPCR) was performed on 39 of the paired tumor and normal lung tissues to validate the mRNA expression of the candidate genes; the samples in the other 6 pairs were insufficient for this analysis. The analysis was carried out on an [Applied Biosystems 7300 Real-Time PCR System](#) (Applied Biosystems, Inc., Foster City, CA, USA) by using Quantinova SYBR Green PCR Master Mix (Qiagen). The relative expression values of tumor and normal lung tissues were calculated by the comparative Ct method. Primers for RT-qPCR are summarized in S4 Table. The expression of β -actin was used as a reference to normalize the other genes' expression. Each experiment was performed in duplicate.

You can see the attachments for the specific experimental process.

 [1076590_hb_quantinova sybr green pcr_1113_ww.pdf](#)

 [rneasy-mini.pdf](#)

 [transcriptor first strand cdna04379012001_en_06.pdf](#)

Western blotting

- 6 Total protein was isolated from lung tissues by lysis of radio immunoprecipitation assay (RIPA) and protected by phenylmethylsulfonyl fluoride (PMSF) following manufacturer's instructions (Beyotime, Shanghai, China). Equal amounts of protein from each sample were separated by 12% sodium dodecyl sulfate–polyacrylamide gel electrophoresis (SDS-PAGE), transferred to polyvinylidene fluoride membranes (Merck Millipore, Billerica, MA, USA) by the semi-dry transfer protocol (Bio-Rad, Hercules, CA, USA), and blocked with 5% skimmed milk. The membranes were then probed with monoclonal antibodies against candidate genes and GAPDH overnight at 4°C. Antibody information is listed in S5 Table. The membranes were then rinsed with TBST and incubated with IgG conjugated to horseradish peroxidase at room temperature for 2 h. The protein bands were visualized on an ImageQuant LAS 500 (GE Healthcare, Pittsburgh, PA, USA) after applying electrochemiluminescent detection reagent (Millipore, Bedford, MA, USA). The intensities of the bands were quantified by using image J software. The band intensity for each candidate gene was subtracted from its own background intensity, normalized against the corresponding GAPDH intensity, and then compared with that of the paired sample. Each experiment was performed in duplicate.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited