

Mixture for M0492 Q5 2X Master Mix

New England Biolabs

Abstract

Citation: New England Biolabs Mixture for M0492 Q5 2X Master Mix. [protocols.io](https://doi.org/10.17504/protocols.io.cijucm)

[dx.doi.org/10.17504/protocols.io.cijucm](https://doi.org/10.17504/protocols.io.cijucm)

Published: 16 Oct 2014

Protocol

Step 1.

Q5 High-Fidelity 2X Master Mix

Step 2.

10 μ M Forward Primer

■ ANNOTATIONS

Florence Servais 28 Jul 2015

Primers from 100 μ M stock => 1/10 dilution first

rPAP1 into secNLuc reporter plasmid: rPAP1secNLucGibsonFor

rPAP1 into CLuc reporter plasmid: rPAP1CLucGibsonFor

STAT3 binding sites + CMV into secNLuc report plasmid: STAT3allsitessecNLucFor

Step 3.

10 μ M Reverse Primer

■ ANNOTATIONS

Florence Servais 28 Jul 2015

From 100 μ M stock => 1/10 dilution first

rPAP1 into secNLuc reporter plasmid: rPAP1secNLucGibsonRev

rPAP1 into CLuc reporter plasmid: rPAP1CLucGibsonRev

STAT3 binding sites + CMV into secNLuc report plasmid: STAT3allsitessecNLucRev

Step 4.

Template DNA

■ ANNOTATIONS

Florence Servais 28 Jul 2015

Use of 0.25ng of each DNA template (plasmids):

#582 (pXP2d2-rPAP1): cc = 0.115 μ g/ μ L => do 1/100 dilution and take **0.217 μ L**

#630 (pSTAT3-CLuc Reporter-MH1): cc = 0.4 μ g/ μ L => do 1/200 dilution and take **0.125 μ L**

Step 5.

Nuclease-Free Water