

Liquid Amplification

Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen

Abstract

For use in "[Obtaining pure cyanophage stock: plaque purification](#)"

Citation: Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen Liquid Amplification. **protocols.io**

dx.doi.org/10.17504/protocols.io.dq65zd

Published: 08 Feb 2016

Guidelines

Use either the liquid amplification method (A) or the plate amplification method (B) in the "[Obtaining pure cyanophage stock: plaque purification](#)" protocol.

Protocol

Step 1.

Add some of the lysate to target host in liquid culture.

Step 2.

After the culture has lysed, remove cell debris via centrifugation.

Step 3.

Filter sterilize the stock.

Step 4.

Store at 4°C until further analysis.