

Elisa  
[Alex Andre](#)<sup>1</sup>  
<sup>1</sup>paris 6  
[dx.doi.org/10.17504/protocols.io.3iagkae](https://doi.org/10.17504/protocols.io.3iagkae)

Working  [Alex Andre](#) 

ABSTRACT  
Titration of antibodies in the serum of patients.

GUIDELINES  
prepare the plate the day before

MATERIALS		
NAME ▾	CATALOG # ▾	VENDOR ▾
PBS		
AggreWell™ 400 24-well plate, 1 pack 1 Plate	34411	Stemcell Technologies
anti-G3BP1 and anti-PABP1 primary antibodies		Abcam

MATERIALS TEXT  
test

SAFETY WARNINGS  
inactivate the serum



BEFORE STARTING  
put the reagents at RT

coating



1  
add ab at [M]1 Micromolar (μM)

2 wash and keep at 4 4 °C

### Blocking

- 3 add 100ul BSA
- 4 incubate  01:00:00  
 20 °C

### Adding of sample

- 5 add the serum in the plate
- 6 incubate  20 °C  
 01:00:00
- 7 Wash

### Detection

- 8 Add Ab at  10 Micromolar (μM)
- 9 incubate  Room temperature  00:00:00



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited