



Oct 09, 2019

Prototype Protocol

George Moberly¹¹GENiE Life Sciences*In Development*[dx.doi.org/10.17504/protocols.io.73yhqpw](https://doi.org/10.17504/protocols.io.73yhqpw)

George Moberly ⚡

ABSTRACT

Prototype Protocol

- 1 Add 200 μ l SWS to the tube. Pipette to mix.
- 2 Place on the 50°C heat block for 30 minutes.
- 3 Immediately place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 4 Remove and discard all supernatant from the tube.
- 5 Remove from the magnetic stand. 1. Add 200 μ l SWS to the tube. Pipette to mix.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited