

# Working Bead Solution

Matthew Sullivan

## Abstract

For use in [Wet-mount Method for Enumeration of Aquatic Viruses](#)

**Citation:** Matthew Sullivan Working Bead Solution. **protocols.io**

[dx.doi.org/10.17504/protocols.io.c84zyv](https://dx.doi.org/10.17504/protocols.io.c84zyv)

**Published:** 07 Jan 2016

## Materials

 Microsphere silica beads (2.34um diameter) [SS04N/4186](#) by Contributed by users

## Protocol

### Step 1.

Dilute stock bead solution 10-fold into PBS to obtain a concentration of 108 beads mL<sup>-1</sup> ; store at 4°C.