

Antibiotic sensitivity testing - *Pyramimonas*, *Bigelowiella*, *Eutreptiella*, *Neovahlkampfia*, *Pseudonitzschia*

Anna Vanclová and Naoji Yubuki

Abstract

Organisms used for DNA transformation experiments were tested for sensitivity to five antibiotics in two different concentrations. Cell culture growing in F2 media with 3% seawater was divided to fifteen 15ml tubes, i.e. three tubes for testing each of the antibiotics: antibiotic was added to two samples to the final concentrations noted in Table 1, the third tube in the set served as a control. The cultures were incubated under 12h light cycle and 20°C and checked by light microscopy 1, 4, 7 and 19 days after the addition of antibiotics.

antibiotic	puromycin	blastidicin	geneticin	hygromycin	zeocin
concentration 1	5 µg/ml	5 µg/ml	200 µg/ml	200 µg/ml	100 µg/ml
concentration 2	10 µg/ml	10 µg/ml	400 µg/ml	400 µg/ml	200 µg/ml

Table 1: Tested antibiotics and concentrations

Citation: Anna Vanclová and Naoji Yubuki Antibiotic sensitivity testing - *Pyramimonas*, *Bigelowiella*, *Eutreptiella*, *Neovahlkampfia*, *Pseudonitzschia*. **protocols.io**
dx.doi.org/10.17504/protocols.io.hh3b38n

Published: 06 Jun 2017

Protocol