May 02,

Transfection of primary endothelial cells with siRNA 👄

PLOS One

2019

Working

Thomas Mudersbach¹, Daniel Siuda¹, Karin Kohlstedt¹, Ingrid Fleming¹

¹Institute for Vascular Signalling, Centre for Molecular Medicine, Goethe University, Frankfurt am Main, Germany

dx.doi.org/10.17504/protocols.io.xtgfnmw



Daniel Siuda

EXTERNAL LINK

https://doi.org/10.1371/journal.pone.0216218

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Epigenetic control of the angiotensin-converting enzyme in endothelial cells during inflammation PLOS ONE

Preparing the endothelial cells for transfection (12-well plate)

Cells were seeded the day before transfection at a density of 3×10^4 cells/cm².

Transfection

Cells were washed twice with PBS and than starved for at least ⊚ 02:00:00 in ⊒500 µI serum-free basal MCDB131 Medium at

8 37 °C

© 00:20:00 before incubation is finshed prepare the transfection mixture (volumes are for a single well of a 12-well plate)

- in one tube (A) add

 1.5 μl of siRNA (50 pmol stock) to

 125 μl Opti-MEM
- in a second tube (B) add 3.75 µl of Lipofectamin RNAiMax into 125 µl Opti-MEM

Mix each tube gently and incubate at room temperatur for © 00:10:00

After incubation the siRNA solution from tube A was added to the lipofectamin solution in tube B and incubated for \odot 00:10:00 at room temperature

■225 µl of the siRNA/lipofectamin solution was added to each well

After © 04:00:00 to © 06:00:00 the medium was changed to endothelial cell growth medium

The endothelial cells were incubated for (§ 48:00:00 at § 37 °C the washed with PBS and snap frozen in liquid nitrogen for further analysis

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited