

BIOL 354W - Research Methods in Advance Microbiology

Version 7

Rosa Leon

Abstract

This protocol series will guide students through the experience of analyzing metagenomic data.

Citation: Rosa Leon BIOL 354W - Research Methods in Advance Microbiology. **protocols.io**

dx.doi.org/10.17504/protocols.io.m26c8he

Published: 05 Feb 2018

Protocol

Introduction to BIOL 354W, sequencing data and bioinformatics

Step 1.

[BIOL 354W Jan 16th](#)

[BIOL 354W Jan 18th](#)

Command line tutorial

Step 2.

In order to do bioinformatics, we first need to become comfortable using the computational language and basic skills that will allow you to analyze data.

Open this link in Chrome

 LINK:

http://rik.smith-unna.com/command_line_bootcamp/

 NOTES

Marcia Smith 29 Jan 2018

change to:

In order to do bioinformatics, we first need to become comfortable using the computational language and basic skills that will allow you to analyze data.

DNA quality assessment and assurance

Step 3.

The first step in analyzing the sequencing data set is to assess the quality of the sequence, and then to edit the dataset in order to retain only the highest quality sequences for the following analysis.

To this end we will use: FastQC - A high throughput sequence QC analysis tool

Familiarize yourself with the software by looking at their [web page](#) - check out the video tutorial!

cmd COMMAND

```
scp -r username@bio-server-2.willamette.edu:/home/username/folder_with_fastqc_file ~/Desktop/
```

Now that the software has run and you have folders and files with data, you should look at the data to assess the quality and make a decision about the quality control step that we will work on next. For this you can unzip your folder where there will be detail information about the results, as well as a summary of the run. You can also download the .html file to look at the graphic representation of the run, the same format you experienced on the FastQC web and tutorial

📌 NOTES

Rosa Leon 14 Jan 2018

You can perform the FastQC file on .fastq files and also in .fastq.gz files or compressed files

Rosa Leon 30 Jan 2018

This step must be done from a Terminal window that is looking at your own computer and not connected to the server

Marcia Smith 29 Jan 2018

Change to:

The first step in analyzing the sequencing data set is to assess the quality of the sequence, and then to edit the dataset in order to retain only the highest quality sequences for the following analysis.

Assuring DNA sequencing quality using Trimmomatic

Step 4.

Trimmomatic: A flexible read trimming tool for Illumina NGS data ([Website](#))

Description

Trimmomatic performs a variety of useful trimming tasks for Illumina paired-end and single ended

data. The selection of trimming steps and their associated parameters are supplied on the command line.

The current trimming steps are:

- ILLUMINACLIP: Cut adapter and other illumina-specific sequences from the read.
- SLIDINGWINDOW: Perform a sliding window trimming, cutting once the average quality within the window falls below a threshold.
- LEADING: Cut bases off the start of a read, if below a threshold quality
- TRAILING: Cut bases off the end of a read, if below a threshold quality
- CROP: Cut the read to a specified length
- HEADCROP: Cut the specified number of bases from the start of the read
- MINLEN: Drop the read if it is below a specified length
- TOPHRED33: Convert quality scores to Phred-33
- TOPHRED64: Convert quality scores to Phred-64

cmd **COMMAND**

```
java -jar /opt/BioInfo_tools/Trimmomatic-0.36/trimmomatic-0.36.jar PE -threads 5 -
phred33 input_forward.fq.gz input_reverse.fq.gz output_forward_paired.fq.gz output_forward_
unpaired.fq.gz output_reverse_paired.fq.gz output_reverse_unpaired.fq.gz ILLUMINACLIP:/opt/
BioInfo_tools/Trimmomatic-0.36/adapters/TruSeq3-
PE.fa:2:30:10 LEADING:15 TRAILING:15 SLIDINGWINDOW:4:15 MINLEN:36
input_forward.fq.gz = " the exact name of your forward or R1 sequence file" input_reverse.fq.gz =
" the exact name of your forward or R2 sequence file" output_forward_paired.fq.gz = write in what
you would like the output to be called Eg. 3A_trimmed_R1_paired.fastq.gz"
output_forward_unpaired.fq.gz = write in what you would like the output to be called Eg.
3A_trimmed_R1_unpaired.fastq.gz" output_reverse_paired.fq.gz = write in what you would like the
output to be called Eg. 3A_trimmed_R2_paired.fastq.gz" output_reverse_unpaired.fq.gz = write in
what you would like the output to be called Eg. 3A_trimmed_R2_unpaired.fastq.gz" Try to run this
command as it is with quality of Q15 (SLIDINGWINDOW:4:15) as currently stated in the command
and then with Q30 (SLIDINGWINDOW:4:30). Record the number % of out put sequences per each.
```

Metagenomic assembly

Step 5.

To assemble our metagenomes we will try two different assemblies and compare them. We will try IDBA_UD and Megahit assemblies. These are going to be one of the most time intensive processes that we will do in the class.

Megahit github - <https://github.com/voutcn/megahit/>

Megahit article - <https://academic.oup.com/bioinformatics/article/31/10/1674/177884>

IDBA_UD - <https://github.com/loneknightpy/idba>

IDBA_UD article - <https://academic.oup.com/bioinformatics/article/28/11/1420/266973>

cmd **COMMAND**

```
/opt/BioInfo_tools/idba/bin/idba_ud -r merged_reads.fa -o output_dir --num_threads 5
```

Once the read files are converted into fasta and in consecutive order then the assembly can be run
merged_reads.fa = your new generated merged fasta sequences files exactly as you called them
output_dir = a folder to store the assembly output, you choose the folder name

Assessing the quality of the assemblies

Step 6.

We can investigate assembly statistics to compare which assembly is best between the two assemblies utilized. For this we can use a software called Quast.

Metrics based only on contigs:

- Number of large contigs (i.e., longer than 500 bp) and total length of them.
- Length of the largest contig.
- N50 (length of a contig, such that all the contigs of at least the same length together cover at least 50% of the assembly).
- Number of predicted genes, discovered either by GeneMark.hmm (for prokaryotes), GeneMark-ES or GlimmerHMM (for eukaryotes), or MetaGeneMark (for metagenomes).

cmd **COMMAND**

```
/opt/BioInfo_tools/quast/metaquast.py contig.fa --gene-finding -t 5
```

QUAST evaluates genome assemblies by computing various metrics.

Binning assembled metagenomes

Step 7.

MaxBin is a software for binning assembled metagenomic sequences based on an Expectation-Maximization algorithm.

Users provide the assembled metagenomic sequences and the reads coverage information or sequencing reads. MaxBin will report genome-related statistics, including estimated completeness, GC content and genome size in the binning summary page.

MaxBin article - <https://academic.oup.com/bioinformatics/article/32/4/605/1744462>

cmd **COMMAND**

```
perl /opt/BioInfo_tools/MaxBin-2.2.4/run_MaxBin.pl -contig "assembled.fa" -  
reads "concatenated reads fasta" -out "out directory"
```

MaxBin requires the assembled contains file and also the file that contains the sequence reads
assembled.fa = the path to your assembled contigs file concatenated reads fasta = the path to
your reads, these reads most all be in one file and concatenate (or paired R1 followed by R2

reads). This you can get from your IDBA fq2fa run out directory = a directory that you create to save your bins

Assessing the quality of your bins via CheckM

Step 8.

Checkm article - <http://genome.cshlp.org/content/25/7/1043>

Also check out the website for information on CheckM - [CheckM website](#)

cmd **COMMAND**

```
/usr/bin/checkm lineage_wf ./bins_folder ./checkm_out_folder
```

CheckM will assess the quality of each of your bins. All bins must be in the same directory/folder.

All bins must have a .fasta ending bins_folder = the path to the folder where your bins are located

Use VizBin to further curate your bins

Step 9.

Perform taxonomic identification using Phylosift

Step 10.