

# **Design of genome-wide HD-FISH probes**

Magda Bienko, Nicola Crosetto, Leonid Teytelman, Sandy Klemm, Shalev Itzkovitz, and Alexander van Oudenaarden

# **Abstract**

This protocol describes the design of primer pairs against the human genome for the synthesis of probes for high-definition DNA FISH (HD-FISH). This pipeline selects PCR primer pairs with optimal thermodynamic features, delimiting amplicons 200–220 nucleotides in length, and filters out primer pairs amplifying multiple targets as well as cross-hybridizing amplicons. Using such primers, highly specific double-stranded probes can be rapidly generated for virtually any desired genomic locus by fluorescently labeling pooled amplicons after PCR.

While this protocol describes the design against the human genome, we have also used it to geneate a genome-wide library for mouse. The design and method should work across other organisms as well.

For more information please see the full paper (and the dedicated <a href="https://hdf.eu/hdfish

**Citation:** Magda Bienko, Nicola Crosetto, Leonid Teytelman, Sandy Klemm, Shalev Itzkovitz, and Alexander van Oudenaarden Design of genome-wide HD-FISH probes. **protocols.io** 

dx.doi.org/10.17504/protocols.io.d6m9c5

Published: 09 Nov 2015

# **Guidelines**

#### **OVERVIEW**

- 1. Each chromosome is split into 500-bp tiled fragments (sliding in 100bp steps)
- 2. Each fragment is searched against the human genome with blat\*
- 3. Unique fragments are merged if overlapping
- 4. Primer 3 is run on the merged sequences for tiled primer pairs design.
- 5. Each primer pair is checked for uniqueness by e-PCR.
- 6. Check each probe for uniqueness against the genome with a more sensitive BLAT search.\*\*

\*First round of BLAT (on 500bp-windows) uses 80% cutoff for matches.

\*\*Final round of BLAT (on the 200bp probes) uses 70% cutoff.

#### **SOFTWARE**

Homology searches are done with BLAT (<a href="http://genome.ucsc.edu/FAQ/FAQblat.html">http://genome.ucsc.edu/FAQ/FAQblat.html</a>)

Primers are designed with "Primer3" (<a href="http://primer3.sourceforge.net/">http://primer3.sourceforge.net/</a>)

# Total sequence for primer design, with single-hits only:

```
1,228,290,200bp (step 3)
5,120,725 initial primers designed (step 4)
5,080,020 after in silico PCR (step 5)
4,823,784 final probes
-Final round of BLAT (on the 200bp probes) uses 70% cutoff.
```

- e-PCR search done with maximum of 2bp mismatch in the primers.

# **NOTES**

The pipeline uses a cluster, so many of the commands are specific to breaking up the analysis files and submitting parallel jobs to the cluster nodes.

All scripts are available at hdfish.edu.

#### **ANALYSIS**

Want to get a feel for the coverage of chromosomes.

- How many probes per 10KB?
- How many deserts with very few probes?
- Excluding the deserts, what is the probe density?

Below are the SQL analyses to answer the above.

LOAD DATA LOCAL INFILE '/home/lenny/projects/hd\_fish/db/round2\_unique\_probes.tbl' INTO TABLE probes;

```
-Count "N"s in the fasta file
(on rous)
cd /home/nylenny/hdfish/databases/fasta/
```

```
for file in chr*.fa; do cat $file | perl -ne 'chomp;if (m/^>/){$chrom=$_;$chrom=\sims/>//;print "$chrom\t";}else{$n_count = tr/nN//;$totaln+=$n_count}END{print "$totaln\n";}' >>
```

```
human_chr_ncounts.txt; done

CREATE TABLE chroms(
   chrom VARCHAR(20),
   size INT
);

CREATE TABLE ncounts(
   chrom VARCHAR(20),
   totalN INT
);
```

LOAD DATA LOCAL INFILE '/home/lenny/projects/hd\_fish/db/human\_chr\_ncounts.tbl' INTO TABLE ncounts;

create table chroms\_full select chroms.chrom,chroms.size,ncounts.totalN n\_count, size-totalN adjusted size from chroms, ncounts where chroms.chrom=ncounts.chrom;

```
drop table chroms;
drop table ncounts;
alter table chroms full rename to chroms;
```

- Count probes/10KB for each chromosome

CREATE TABLE temp (select chrom,count(\*) total probes from probes group by chrom);

create table probes\_bychrom select pb.chrom,pb.total\_probes,chroms.size chrom\_size, round(10000\*pb.total\_probes/chroms.size,1) probes\_per10kb from temp pb, chroms WHERE pb.chrom=chroms.chrom order by chroms.size DESC;

create table probes\_bychrom\_adj select pb.chrom,pb.total\_probes,chroms.adjusted\_size adj\_chrom\_size, round(10000\*pb.total\_probes/chroms.adjusted\_size,1) probes\_per10kb from temp pb, chroms WHERE pb.chrom=chroms.chrom order by chroms.size DESC;

drop table temp;

Look at 10KB-windows
 create table windows select chrom,round(start/10000) start,count(\*) total from probes group by chrom,round(start/10000);

# **Protocol**

#### **BLAT** database setup

# Step 1.

Download the human genome

```
cmd COMMAND
```

cd /home/nylenny/hdfish/databases/fasta

ftp://hgdownload.cse.ucsc.edu/goldenPath/currentGenomes/Homo sapiens/bigZips/chromFa.tar.gz

# BLAT database setup

#### Step 2.

Remove all strange .fa files accept for chr1.fa chr2.fa...

```
cmd COMMAND
rm chrUn_gl0002* *_gl* *hap*
```

#### **BLAT** database setup

#### Step 3.

Make the 2bit file for BLAT

# NOTES

# Lenny Teytelman 09 Nov 2015

"qsub" is a command for submitting jobs to cluster nodes using the <u>PBS Pro</u> software (Portable Batch System).

#### **BLAT** database setup

# Step 4.

Make an 11.ooc file

```
cmd COMMAND
cd /home/nylenny/hdfish/databases/blatdb
ls -1 *.2bit > human_genome_chrom.list
~/programs/blat/blat -makeOoc=11.ooc -
tileSize=11 human_genome_chrom.list /dev/null /dev/null
```

# Split chroms into 500-bp fragments

#### Step 5.

Make a fasta file of 500-bp windows, in 100bp slidingsteps for each chromosome.

```
cmd COMMAND
cd /home/nylenny/hdfish/databases/fasta

for file in *.fa
do
        export infile=/home/nylenny/hdfish/databases/fasta/$file
        export outfile=/home/nylenny/hdfish/databases/split_chroms/${file%.fa}._windows.fa
        qsub -o /home/nylenny/hdfish/cluster_eo/ -e /home/nylenny/hdfish/cluster_eo/ -
V /home/nylenny/hdfish/qsub_runs/split_chroms.sh
done
```

#### Split chroms into 500-bp fragments

#### Step 6.

Get human chromosome sizes

```
cmd COMMAND
for file in *.fa;do echo $file; grep -v '>' $file | wc;done > human_chr_sizes.txt

cat human_chr_sizes.txt | perl -
ne 'chomp;if(m/fa$/){s/.fa//;print "$_\t";}else{s/^\s*//;my($lines,$words,$total_chars)=spl
it(" +");$chrom_size=$total_chars-$lines;print "$chrom_size\n";}' > human_chr_sizes.txt2

mv human_chr_sizes.txt2 human_chr_sizes.txt
• NOTES
```

# Lenny Teytelman 09 Nov 2015

Steps 6-8 are a sanity check. Making sure that the resulting split\_windows chrom sizes match the actual size of the human chromosomes in the initial .fasta files.

#### Split chroms into 500-bp fragments

#### Step 7.

Get the sizes of the split chroms

```
cmd COMMAND
cd /home/nylenny/hdfish/databases/split_chroms

for file in *.fa; do tail -2 $file | head -1; done | tr -d '>' |tr '_' '\t' | cut -f 1,3|sort > split chrom sizes.txt
```

# Split chroms into 500-bp fragments

#### Step 8.

Visually inspect the output from this step to compare real chroms against the split ones

```
cmd COMMAND
join -j 1 ../fasta/human_chr_sizes.txt split_chrom_sizes.txt | tr ' ' '\t' | perl -
ne 'chomp;my($chrom,$real_size,$split_size)=split("\t");print "$chrom\t",$real_size-
$split_size,"\n";'
```

#### BLAT, first round

#### Step 9.

BLAT the split windows for each chromosome against the entire human genome, one chromosome at a time.

# NOTES

#### Lenny Teytelman 09 Nov 2015

"qsub" is a command for submitting jobs to cluster nodes using the <u>PBS Pro</u> software (Portable Batch System).

#### Lenny Teytelman 09 Nov 2015

This first round of BLAT (on 500bp-windows) is more permissive and uses 80% cutoff for matches.

# BLAT, first round

#### Step 10.

Check for completion of the blat against each chromosome. (The last blat hit should be close to the end of the guery file.)

```
cmd COMMAND
cd /home/nylenny/hdfish/blat_results

for file in chr3*; do echo $file |perl -
ne 'chomp;s/chr4._windows.fa_to_//;s/.fa.2bit.psl//;print "$_\t";'; tail -1 $file | cut -
f 10 ;done

for file in chr*; do echo $file |perl -
```

```
ne 'chomp;s/chr5._windows.fa_to_//;s/.fa.2bit.psl//;print "\_\t";'; tail -1 $file | cut -f 10 ;done > sizes.txt
```

# Find unique fragments

#### **Step 11.**

Merge all the .psl files for each chromosome

```
cmd COMMAND
qrsh

cd /home/nylenny/hdfish/blat_results
mkdir by_chrom
mv *.psl by_chrom/

for file in chr1 chr10 chr11 chr12 chr13 chr14 chr15 chr16 chr17 chr18 chr19 chr2 chr20 chr
21 chr22 chr3 chr4 chr5 chr6 chr7 chr8 chr9 chrX chrY
do
    qsub -e /home/nylenny/hdfish/cluster_eo/ -
o /home/nylenny/hdfish/blat_results/$file.psl -
b y cat /home/nylenny/hdfish/blat_results/by_chrom/${file}._*.psl
done
```

for file in chr11 chr12 chr13 chr14 chr15 chr16 chr17 chr18 chr19 chr2 chr20 chr21 chr22 chr3 chr4 chr5 chr6 chr7 chr8

#### NOTES

# Lenny Teytelman 09 Nov 2015

In this section, we are processing the BLAT results to identify all 500-bp fragments that have a single hit to the genome.

# Find unique fragments

#### **Step 12.**

Count total BLAT hits per fragment

```
cmd COMMAND
cd /home/nylenny/hdfish/blat_results/

for file in *.psl
do
        export infile=/home/nylenny/hdfish/blat_results/$file
        export outfile=/home/nylenny/hdfish/blat_results/counts/${file%.psl}.blatcounts.txt
        qsub -o /home/nylenny/hdfish/cluster_eo/ -e /home/nylenny/hdfish/cluster_eo/ -
V /home/nylenny/hdfish/qsub_runs/count_blat_hits.sh
done
```

# Merge unique fragments

#### **Step 13.**

Sort by positional start

```
cmd COMMAND
qrsh
cd /home/nylenny/hdfish/blat_results/counts

for file in *.txt
do
    tr '_' '\t' < $file | cut -f 2,4 | sort -k 1n > ${file%.blatcounts.txt}_sorted_blatcounts.csv
done
```

# NOTES

# Lenny Teytelman 09 Nov 2015

In this section, we merge all overlapping single-hit windows and get their DNA sequences.

# Merge unique fragments

#### **Step 14.**

Make a list of windows with just a single hit

```
cmd COMMAND
for file in *.csv
do
  cat $file | perl -
ne 'chomp;my($start,$count)=split("\t");print "spacer\tspacer\t$start\t",$start+499,"\n" if
  $count
```

## Merge unique fragments

# **Step 15.**

Merge the overlapping windows

# NOTES

#### Lenny Teytelman 09 Nov 2015

[We're now done with the round1 blat .psl output files. Can zip them to save space.]

cd /home/nylenny/hdfish/blat\_results/

```
for file in *.psl do qsub -e /home/nylenny/hdfish/cluster_eo/ -o /home/nylenny/hdfish/cluster_eo/ -b y gzip /home/nylenny/hdfish/blat_results/$file done
```

## Merge unique fragments

#### **Step 16.**

On local machine, make a fasta file with the merged unique sequences

```
cmd COMMAND
cd /media/Elements/projects/hdfish/unique_regions

for file in *.txt; do perl ~/projects/hd_fish/scripts/name_intervals.pl -
i $file > ${file%_singles_merged.txt}_unique_toget.txt;done

for file in *toget.txt; do chrom=${file%_unique_toget.txt}; perl ~/programs/Scripts/get_fa
sta_sequence.pl -f /media/Elements/projects/hdfish/databases/fastas/human/${chrom}.fa -
seq_list $file > /media/Elements/projects/hdfish/databases/fastas/unique_regions/${chrom}_u
nique.fa; done
```

#### Design tiled primer pairs for PCR

#### **Step 17.**

Convert the fasta sequences to primer3 input records

```
cmd COMMAND
cd /media/Elements/projects/hdfish/databases/fastas/unique_regions/
for file in *.fa
do
```

```
perl ~/projects/hd_fish/scripts/make_primer3_records.pl -
i $file > /media/Elements/projects/hdfish/databases/primer3/${file%_unique.fa}_unique_for_
primer3.txt
done
```

#### NOTES

# **Lenny Teytelman** 09 Nov 2015

In this section, we run Primer3 to design tiled PCR primer pairs against the unique genome sections identified above.

# Design tiled primer pairs for PCR

# **Step 18.**

Transfer all the files to the cluster: rous:/home/nylenny/hdfish/databases/primer3

# Design tiled primer pairs for PCR

# Step 19.

Run Primer3 on each chrom\_unique file

## Select unique primer pairs

#### Step 20.

Process the Primer3 output.

```
cmd COMMAND

cd /home/nylenny/hdfish/primers/primer3_output

for file in *.txt

do
    perl /home/nylenny/hdfish/scripts/parse_primer3.pl -i $file > $file.parsed

done
```

#### NOTES

#### **Lenny Teytelman** 09 Nov 2015

In this section, we use e-PCR to ensure that each primer pair has a unique match in the genome and will not cross-hybridize.

# Select unique primer pairs

#### Step 21.

Count primers per chromosome (normalized to 10kb)

```
cmd COMMAND
wc *.parsed | tr -s ' ' '\t' | grep -v total | cut -f 2,5 | cut -d '_' -f 1 | tr -
d 'chr' | sort -k 2n > primer_pair_counts.tmp

cat /home/nylenny/hdfish/databases/fasta/human_chr_sizes.txt | grep -v chrM | tr -
d 'chr' | sort -k 1n > human_sizes.tmp

join -1 2 -2 1 primer_pair_counts.tmp human_sizes.tmp | tr -s ' ' '\t' | perl -
ne 'chomp;my($chrom,$primers,$size)=split("\t");print "$chrom\t$primers\t$size\t",int(10000)
```

## Select unique primer pairs

#### Step 22.

```
Make .sts files for e-PCR
  cmd COMMAND
  cd /home/nylenny/hdfish/primers/primer3_output
  for file in *.parsed
  do
   cat $file | perl -
  ne 'chomp;my($id,$leftp,$rightp,$left_start,$right_start)=split("\t");my($window,$wstart,$w
  stop)=split("_",$id);print "$id\t$leftp\t$rightp\t",$wstart+$left_start,"\t",$wstart+$right
  _start,"\n";' | sort -
  k 4n > /home/nylenny/hdfish/primers/epcr_input/${file%_primer3_output.txt.parsed}_primer3_
  rawcoord.txt
  done
  cd /home/nylenny/hdfish/primers/epcr input/
  for file in *.txt
  dΛ
    cat $file | perl -ne 'chomp; my($id,$
Select unique primer pairs
```

# Step 23.

Run e-PCR against the entire human genome

```
cmd COMMAND
cd /home/nylenny/hdfish/primers/epcr_input
export stsfile=chr22 primer3 rawcoord.sts
cd /home/nylenny/hdfish/databases/fasta
for chromosome in *.fa
do
    export chrom=$chromosome
    qsub -V -
o /home/nylenny/hdfish/primers/epcr output/${stsfile% primer3 rawcoord.sts} to $chrom.epcr
e /home/nylenny/hdfish/cluster_eo/${stsfile}_to_$chrom.e /home/nylenny/hdfish/qsub_runs/epc
r.sh
done
```

#### Select unique primer pairs

# Step 24.

Check for completion of batch jobs for each chromosome

```
cmd COMMAND
cd /home/nylenny/hdfish/primers/epcr_output
for file in chr22_*; do echo $file |perl -
ne 'chomp;s/chr3_to_//s;s/.fa.epcr//;print "$_\t";'; tail -1 $file ;done | tr -
s ' ' '\t' | cut -f 1,4
qrsh
cd /home/nylenny/hdfish/primers/epcr output/
```

```
mkdir bychrom
mv * bychrom/
```

for file in chr1 chr10 chr11 chr12 chr13 chr14 chr15 chr16 chr17 chr18 chr19 chr2 chr20 chr 21 chr22 chr3 chr4 chr5 chr6 chr7 chr8 chr9 chrX chrY

qsub -e /home/nylenny/hdfish/cluster\_eo/ -o /home/nylenny/hdfish/primers/epcr\_output/\$fi

# Select unique primer pairs

# Step 25.

Exclude primers where e-PCR match is not to the exact same position as targeted

```
cmd COMMAND
cat chr*.txt| tr -s ' ' '\t' | tr -s '.' '\t' | tr -s ',' '\t' | tr -s '_' '\t' | perl -
ne 'chomp;$line=$_;my($echrom,$estart,$estop,$chrom,undef,undef,$start,$stop)=split("\t",$l
ine);if ($echrom ne $chrom or $start!=$estart or $stop!=$estop){print "$chrom\t$start\t$sto
p\n";}' |uniq|sort|uniq > wrong_match.txt

perl /home/nylenny/hdfish/scripts/exclude_primers_bylist.pl -
i ../epcr_input/all_primers.txt -e wrong_match.txt > unique_primers.txt
```

# BLAT, round 2

# Step 26.

Now screen the probe sequence of the unique primers pairs by BLAT for uniqueness in the genome.

#### NOTES

#### **Lenny Teytelman** 09 Nov 2015

This second round of BLAT (on the 200bp probes) is stricter and uses 70% cutoff to avoid cross-hybridizing probes.

# BLAT, round 2

# Step 27.

Count hits per probe and select single-hit ones only

```
cmd COMMAND
cd hdfish/blat_results/round2/
cat probes_round2_to_chr*.psl > round2_to_allhuman.psl
perl /home/nylenny/hdfish/scripts/count_blat_hits.pl -
i round2_to_allhuman.psl > counts/round2_all_blatcounts.txt
cat round2_all_blatcounts.txt | perl -
```

ne 'chomp;my( $n=1;' > round2\_uniq ue\_probes.txt$ 

**ANNOTATIONS** 

Lenny Teytelman 10 Apr 2017