



May 09,  
2019

Working

## U Mass - Acute lipid infusion [↗](#)

Jason Kim<sup>1</sup>

<sup>1</sup>University of Massachusetts

dx.doi.org/10.17504/protocols.io.xm5fk86

Mouse Metabolic Phenotyping Centers

Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)

Lili Liang

### ABSTRACT

#### Summary:

Triglyceride emulsion and heparin will be intravenously infused for 5 hours to acutely raise circulating fatty acids levels in awake mice. Acute lipid infusion is shown to cause insulin resistance in peripheral organs.

### EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=153&docType=Protocol>

### MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
INTRALIPID® 20%	NDA18-449/S-039, NDA17-643/S-072	Baxter Healthcare		
Heparin	NDC0409-2723-01	Hospira(Pfizer)		
Glycerol	G5516	Sigma Aldrich		

### MATERIALS TEXT

#### Note:

Baxter [RRID:SCR\\_003974](#)

Sigma-Aldrich [RRID:SCR\\_008988](#)

Hospira [RRID:SCR\\_003985](#)

- 1 Survival surgery is performed to establish a chronic indwelling catheter at 5~6 days prior to experiment for intravenous infusion. (refer to M1023: Surgery-jugular vein cannulation)
- 2 Mice are fasted overnight (~15 hours) or for 5 hours prior to the start of experiment.
- 3 Place a mouse in a rat-size restrainer with its tail tape-tethered at one end.
- 4 Expose and flush the intravenous catheter using saline solution. Then, connect the catheter to the CMA Microdialysis infusion pump.
- 5 Collect plasma sample (10 µl) before the start of infusion (basal-0 min) to measure basal fatty acids levels.

- 6 Start the experiment by turning on the pump and intravenously infusing 20% Intralipid (triglyceride emulsion) at 2.5 ml/kg/hr and heparin at 6 U/hr in awake mice.
- 7 For control experiment, 20% glycerol is infused at 2.5 ml/kg/hr in awake mice.
- 8 Collect plasma samples (10 µl each) at 120, 240, and 300 min to measure serum FFA levels.
- 9 At the end of experiment, tissues may be collected for biochemical and molecular analyses.
- 10 Alternatively, a 2-hr hyperinsulinemic-euglycemic clamp may follow this 5-hr acute lipid infusion experiment to measure the effects of elevated fatty acids on insulin sensitivity.  
(refer to M1001: Hyperinsulinemic-euglycemic clamp)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited