

Q5 PCR Mixture for M0302

New England Biolabs

Abstract

Citation: New England Biolabs Q5 PCR Mixture for M0302. [protocols.io](https://doi.org/10.17504/protocols.io.cqqvzv)

[dx.doi.org/10.17504/protocols.io.cqqvzv](https://doi.org/10.17504/protocols.io.cqqvzv)

Published: 29 Jan 2015

Protocol


Step 1.

Q5® Hot Start High-Fidelity 2X Master Mix (M0494S)

 [AMOUNT](#)

25 µl Additional info:

 [REAGENTS](#)

 Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns [M0494S](#) by [New England Biolabs](#)

Step 2.

10 µM Forward Primer, **2.5 µl**

 [AMOUNT](#)

3 µl Additional info:

Step 3.

10 µM Reverse Primer, 2.5 µl

 [AMOUNT](#)

3 µl Additional info:

Step 4.

Template DNA

Step 5.

Nuclease-free water to **50 µl**