NEB Immunoprecip. Buffer

New England Biolabs

Abstract

Citation: New England Biolabs NEB Immunoprecip. Buffer. protocols.io

dx.doi.org/10.17504/protocols.io.crbv2m

Published: 02 Feb 2015

Protocol

Step 1.

50 mM NaCl

Step 2.

10 mM Tris-HCl (pH 7.4)

Step 3.

1 mM EDTA

Step 4.

1 mM EGTA (pH 8.0)

Step 5.

0.2 mM sodium ortho-vanadate

Step 6.

0.2 mM PMSF

Step 7.

1% Triton X-100

Step 8.

0.5% NP-40