



Demultiplexing Nanopore reads with LAST [↗](#)

David Eccles¹

¹Malaghan Institute of Medical Research (NZ)

[dx.doi.org/10.17504/protocols.io.xj3fkqn](https://doi.org/10.17504/protocols.io.xj3fkqn)

MinION user group for high molecular weight DNA extraction from all kingdoms



David Eccles

Malaghan Institute of Medical ...



ABSTRACT

This protocol is for a semi-manual method for read demultiplexing, as used after my presentation *Sequencing DNA with Linux Cores and Nanopores* to work out the number of reads captured by different barcodes.

Input: reads as a FASTQ file, barcode sequences as a FASTA file

Output: reads split into single FASTQ files per target [barcode]

Note: barcode / adapter sequences are not trimmed by this protocol

EXTERNAL LINK

<https://doi.org/10.5281/zenodo.2535894>

PROTOCOL STATUS

In development

We are still developing and optimizing this protocol

Generating Barcode Index

- 1 Prepare a FASTA file containing barcode sequences (see attached FASTA file). To reduce the chance of mismatched adapters, this should *only* contain the barcode sequences. That restriction means this approach will not work for short reads, where the barcode sequences are very likely to occur within sequences.

☐ [barcode_base.fa](#)

- 2 Prepare the LAST index for the barcode file. This will generate seven additional files of the form <index name>.XXX:

```
lastdb barcode_base.fa barcode_base.fa
```

Mapping Reads to Barcodes

- 3 Combine all input reads into a single file

```
pv ../called_all/*.fastq | gzip > reads_all.fastq.gz
```

Note: I'm using the pipe viewer command *pv* to produce a progress indicator while the command is running. If this command is not available, it can be replaced with *cat* with no change in function (apart from not showing progress).

Use LAST in FASTQ alignment mode (-Q 1) to map the reads. In this example, it is distributed over 10 processing threads (-P 10). Here *maf*-

- 4 *convert* is used to convert to a single line per match, *cut* retains only the barcode and read IDs, and *uniq* is used to make sure that multiple same barcodes per read (e.g. for reverse / complement barcodes at each end) will not produce duplicates:

```
lastal -Q 1 -P10 barcode_base.fa <(pv reads_all.fastq.gz) | \
maf-convert tab | grep -v '^#' | cut -f 2,7 | uniq | \
gzip > barcode_assignments.txt.gz
```

For a more stringent search, the output of *lastal* can be piped through *last-map-probs*, which will reduce the likelihood of a partial barcode match to other DNA sequences. The downside is that this is more likely to drop reads due to slight mismatches in the barcode portion of the read:

```
lastal -Q 1 -P10 barcode_base.fa <(pv reads_all.fastq.gz) | last-map-probs | \
maf-convert tab | grep -v '^#' | cut -f 2,7 | uniq | \
gzip > barcode_assignments.txt.gz
```

The output of this command will be a gzipped tab-separated 2-column file with barcode names in the first column, and read IDs in the second column.

Optional [but recommended]: filtering chimeric reads

- 5 Identify reads with multiple barcodes (i.e. potentially chimeric reads). The *sort* command sorts by the second field (read ID), then *uniq* identifies duplicated lines when ignoring the first field (barcode).

```
pv barcode_assignments.txt.gz | zcat | sort -k 2,2 | \
uniq -f 1 -D | gzip > duplicate_assignments.txt.gz
```

- 6 Exclude duplicate read IDs from the fastq file. This uses one of my own scripts, [fastx-fetch.pl](#), to do this directly from a FASTQ file.

```
~/scripts/fastx-fetch.pl -v -i <(zcat duplicate_assignments.txt.gz | \
awk '{print $2}') <(pv reads_all.fastq) | gzip > reads_noChimeric.fastq.gz
```

Splitting Read File Per Barcode

- 7 Create a file containing barcode read counts

```
pv barcode_assignments.txt.gz | zcat | awk '{print $1}' | \
sort | uniq -c > barcode_counts.txt
```

- 8 For each discovered barcode, find the corresponding read IDs, then extract those IDs out of the read FASTQ file. This uses one of my own scripts, [fastx-fetch.pl](#), to do this directly from a FASTQ file:

```
for bc in $(awk '{print $2}' barcode_counts.txt);
do echo "*** ${bc} ***";
fastx-fetch.pl -i <(zgrep ${bc} barcode_assignments.txt.gz | awk '{print $2}') \
<(pv reads_all.fastq.gz) | \
gzip > reads_${bc}.fastq.gz;
done
```

Or, if working from the non-chimeric reads:

```
for bc in $(awk '{print $2}' barcode_counts.txt);
do echo "*** ${bc} ***";
fastx-fetch.pl -i <(zgrep ${bc} barcode_assignments.txt.gz | awk '{print $2}') \
<(pv reads_noChimeric.fastq.gz) | \
```

```
gzip > reads_${bc}.fastq.gz;  
done
```

Note: this step processes through the read file once per barcode, which could take a while depending on how many barcodes are detected.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited