

M0368 RNA and Primer mixture

New England Biolabs

Abstract

Citation: New England Biolabs M0368 RNA and Primer mixture. **protocols.io**
dx.doi.org/10.17504/protocols.io.chut6v
Published: 07 Oct 2014

Protocol

Step 1.

RNA solution 1ng-1 µg total RNA or 1-100 ng polyA-selected RNA

Step 2.

Primer (50uM dT or 60uM Random Primer Mix) **2 µL**

Step 3.

10mM dNTP mix **1 µL**

Step 4.

Nuclease-free H₂O to **final volume of 12 µL**