

Non-Centrifuged Sample Steps

Matthew Sullivan

Abstract

For [One-step growth curves for Cellulophaga phages protocol](#).

Citation: Matthew Sullivan Non-Centrifuged Sample Steps. **protocols.io**

dx.doi.org/10.17504/protocols.io.ddk24v

Published: 02 Feb 2016

Protocol

Step 1.

Pipet 100 μ l from the flask into 900 μ l of MSM in a 1.5 ml tube

■ ANNOTATIONS

Bonnie Poulos 15 Mar 2016

Refer to the protocol "One-step growth curves for Cellulophaga phages" for details regarding the preparation and dilution of the phage preparation.

Step 2.

Vortex briefly

Step 3.

Plate 100 μ l