



Dec 09, 2019

pL0R-LacZ

Bernardo Pollak¹¹J. Craig Venter Institute/PUC

1

Works for me

dx.doi.org/10.17504/protocols.io.baa8iahw

uLoop - Diatom Parts Registry



Bernardo Pollak

Millennium Institute for Integrative Biology - iBio



ABSTRACT

Name: pL0R-LacZ

Description: a L0 receiver vector which can be used for domestication using SapI-mediated restriction ligation or Gibson assembly. Contains lacZ for blue/white colony screening. Requires X-GAL, IPTG optional depending on E. coli strain.

Level 0

Overhangs: Special (cuts into BsaI sites)

Length: 2646

Antibiotic selection: Spectinomycin

Other info: pMB1 ori

Characterisation: ++++ (substantial), works as expected.

Primers:

U1F	CATTACTCGCATCCATTCTC
U1R	GAGACGAGACGAGACAGCCT
UXF	CCAGGATACATAGATTACCA
UXR	GGTGAAGGGCTCGGAGTTG
L0F	TGGGCTGCCTGTATCGAGTG
L0R	GCTGGCCTTTTGCTCACA

L0 F works great for sequencing.

Submitter: Bernardo Pollak

Cloned by: Bernardo Pollak

Availability: Upon request at the Pollak/Dupont/Federici labs at the moment and will be available soon from Addgene.

EXTERNAL LINK

<https://www.uloop.org/domestication.html>

pL0R-lacZ.gb



This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited