

Stranded Mapping from Long Reads version 2

David Eccles

Abstract

This protocol demonstrates how to convert raw long reads produced using a strand-specific sequencing protocol (e.g. ONT's strand-switching protocol) into strand-specific mapped reads.

The general idea is to use LAST to identify the adapter orientation relative to the genome, and then use that information to split BAM files up and recombine them to create two strand-specific files that are displayable in a genome browser.

Citation: David Eccles Stranded Mapping from Long Reads. **protocols.io**

dx.doi.org/10.17504/protocols.io.n8ddhs6

Published: 04 Apr 2018

Before start

I have written [my own script](#) to process LAST results into a CSV format, which makes it easier to do line-by-line data filtering. I have also created a [fastq filtering script](#) that helps for filtering reads into different files.

You will also need access to the following free and open-source software programs:

- [LAST](#)
- [minimap2](#)
- [samtools](#)

And the following additional data files:

- a FASTA file containing strand-specific primer / adapter sequences.

Protocol

Read Correction

Step 1.

Download canu v1.6 or later

Read Correction

Step 2.

I prefer starting off my data analysis with a read correction with Canu (ideally v1.7 now that it's out, because that attempts correction of *all* reads, but here I use canu v1.6). I use minimap as the mapper to speed this up. The genomeSize parameter should be approximately a tenth to a fortieth of the number of bases in your dataset to make sure that no sequences are excluded (bigger is better, as long as Canu doesn't freak out about memory consumption):

(This creates a file 4T1_BC06/4T1_BC06.correctedReads.fasta.gz.)

cmd COMMAND

```
~/install/canu/canu-1.6/Linux-  
amd64/bin/canu overlapper=minimap genomeSize=100M minReadLength=100 minOverlapLength=30 -  
correct -p 4T1_BC06 -d 4T1_BC06 -nanopore-raw workspace/pass/barcode06/fastq_runid_*.fastq  
Read correction using Canu v1.6
```

Chimeric read filtering

Step 3.

The next step I carry out is a basic read-level QC to exclude [chimeric reads](#). [Porechop](#) can be used for this, although that removes adapters by default, which is not particularly useful in this case.

I use LAST to search for adapter sequences within the corrected reads, pass it through [my conversion script](#), and extract out duplicated mappings (i.e. where the same read/adaptor pair appears more than once in the mapping results):

cmd COMMAND

```
lastal -  
P 10 ONT_barcode_adapters.fa <(zcat 4T1_BC06/4T1_BC06.correctedReads.fasta.gz) | ~/scripts  
/maf_bcsplit.pl | awk -F',' '{print $1,$2}' | sort | uniq -  
d | awk '{print $1}' | uniq > reads_with_duplicated_adapters.txt  
Mapping adapter sequences with LAST
```

Chimeric read filtering

Step 4.

I use this file to filter out chimeric reads from the corrected dataset using another [fastq filtering script](#) I've created:

cmd **COMMAND**

```
~/scripts/fastx-fetch.pl -v -  
i reads_with_duplicated_adapters.txt 4T1_BC06/4T1_BC06.correctedReads.fasta.gz | gzip > 4T1_BC06.correctedReads.uniqueOnly.fasta.gz  
Removing reads with duplicated adapters
```

Adapter filtering

Step 5.

Along roughly similar lines to the chimeric read filtering, I then look for the strand switch and VNT adapters in the sequences. For a forward-orientated query, I expect the strand switch primer to be in the forward direction, and the VNP primer to be in the reverse direction.

Adapter filtering

Step 6.

Transcript-forward reads are identified. In this case a sequence is transcript-forward if the `ONT_SSP` primer is in the forward orientation, or the `ONT_VNP` primer is in the reverse orientation.

cmd **COMMAND**

```
lastal -  
P 10 ONT_barcode_adapters.fa <(zcat 4T1_BC06.correctedReads.uniqueOnly.fasta.gz) | ~/scripts/maf_bcsplit.pl | grep -e 'ONT_SSP,+' -e 'ONT_VNP,-' | awk -  
F',' '{print $1}' | sort | uniq > fwdQry_seqs_BC06.txt  
Identifying transcript-forward reads with LAST
```

Adapter filtering

Step 7.

Transcript-reverse reads are identified. In this case a sequence is transcript-reverse if the `ONT_SSP` primer is in the reverse orientation, or the `ONT_VNP` primer is in the forward orientation.

cmd **COMMAND**

```
lastal -  
P 10 barcodes_primerSeqs.fa <(zcat 4T1_BC06.correctedReads.uniqueOnly.fasta.gz) | ~/scripts/maf_bcsplit.pl | grep -e 'ONT_SSP,-' -e 'ONT_VNP,+' | awk -  
F',' '{print $1}' | sort | uniq > revQry_seqs_BC06.txt  
Identifying transcript-reverse reads with LAST
```

Adapter filtering

Step 8.

The read name lists are used by my [read filtering script](#) to split the original reads up into transcript-forward and transcript-reverse subsets:

cmd **COMMAND**

```
~/scripts/fastx-fetch.pl -  
i fwdQry_seqs_BC06.txt 4T1_BC06.correctedReads.uniqueOnly.fasta.gz | gzip > fwd_4T1_BC06.correctedReads.uniqueOnly.fasta.gz
```

```
~/scripts/fastx-fetch.pl -  
i revQry_seqs_BC06.txt 4T1_BC06.correctedReads.uniqueOnly.fasta.gz | gzip > rev_4T1_BC06.co  
rrectedReads.uniqueOnly.fasta.gz  
Filtering reads into transcript-forward and transcript-reverse subsets
```

Read Mapping

Step 9.

Download minimap2 v2.10 or later

 **SOFTWARE PACKAGE (Linux)**

minimap2, 2.10 

Heng Li
<https://github.com/lh3/minimap2>

Read Mapping

Step 10.

Now that the reads have been oriented, the mapping can be done. I use minimap2 for mapping long reads to a transcriptome.

cmd **COMMAND**

```
~/install/minimap2/minimap2 -t 10 -a mmus_ucsc_all_cdna.idx -  
x splice fwd_4T1_BC06.correctedReads.uniqueOnly.fasta.gz > fwd_4T1_BC06.CU_vs_mmus.bam
```

```
~/install/minimap2/minimap2 -t 10 -a mmus_ucsc_all_cdna.idx -  
x splice rev_4T1_BC06.correctedReads.uniqueOnly.fasta.gz > rev_4T1_BC06.CU_vs_mmus.bam
```

Mapping reads to a transcriptome using minimap2

Splitting and recombining strands

Step 11.

Download samtools v1.8 or later

 **SOFTWARE PACKAGE (Linux)**

SAMtools, 1.8 

Wellcome Trust Sanger Institute
<https://github.com/samtools/samtools>

Splitting and recombining strands

Step 12.

Mapped BAM files are split based on how they are oriented to the *genome* using SAMtools. This uses the ``-F`` or ``-f`` options of ``samtools view`` to exclude or include reads (respectively) that are mapped to the genome in the forward or reverse orientation (respectively).

cmd **COMMAND**

```
samtools view -F 0x10 -b fwd_4T1_BC06.CU_vs_mmus.bam > fwd_fwd_4T1_BC06.CU_vs_mmus.bam
```

```
samtools view -f 0x10 -b fwd_4T1_BC06.CU_vs_mmus.bam > rev_fwd_4T1_BC06.CU_vs_mmus.bam
```

```
samtools view -F 0x10 -b rev_4T1_BC06.CU_vs_mmus.bam > fwd_rev_4T1_BC06.CU_vs_mmus.bam
```

```
samtools view -f 0x10 -b rev_4T1_BC06.CU_vs_mmus.bam > rev_rev_4T1_BC06.CU_vs_mmus.bam
```

Splitting BAM files with samtools based on genome-relative orientation

Splitting and recombining strands

Step 13.

And finally the files are recombined to identify forward and reverse-encoded transcripts via `samtools merge`. I'm calling the "pos" strand the one that is encoded in the same direction as the genome, and the "neg" strand the reverse-complement direction:

cmd **COMMAND**

```
samtools merge pos_4T1_BC06.CU_vs_mmus.bam rev_fwd_4T1_BC06.CU_vs_mmus.bam fwd_rev_4T1_BC06.CU_vs_mmus.bam
```

```
samtools merge neg_4T1_BC06.CU_vs_mmus.bam fwd_fwd_4T1_BC06.CU_vs_mmus.bam rev_rev_4T1_BC06.CU_vs_mmus.bam
```

Merging forward and reverse-oriented reads with ``samtools merge``

Splitting and recombining strands

Step 14.

NOTES

David Eccles 04 Apr 2018


Alternatively, it's possible to set the BAM flags for read 1 and read 2 for the reverse and forward-encoded adapters respectively, which allows everything to then be combined into a single BAM file (and treated in the same way as a strand-specific Illumina BAM file).

Sanity Check

Step 15.

EXPECTED RESULTS

If this has worked properly, then mapping human or mouse to the mitochondrial genome should show most expression appearing on the positive strand, with a small scattering of negative-strand expression, a bit like this:



Warnings

These scripts have been slightly modified from scripts that I have run. Consider them demonstrative: pay attention to the words rather than the script.