

Mixture for M0290 CIP

New England Biolabs

Abstract

Citation: New England Biolabs Mixture for M0290 CIP. [protocols.io](https://doi.org/10.17504/protocols.io.chrt55)

[dx.doi.org/10.17504/protocols.io.chrt55](https://doi.org/10.17504/protocols.io.chrt55)

Published: 09 Oct 2014

Protocol

Step 1.

DNA $\geq 1 \mu\text{l}$

■ ANNOTATIONS

Peter Rüthemann 14 Oct 2015

Cas9

Step 2.

Restriction Enzyme Buffer (10X) $2 \mu\text{l}$

Step 3.

Restriction Endonuclease $1 \mu\text{l}$

Step 4.

H₂O, purified to $20 \mu\text{l}$