

# Mixture for PCR Amplification (E7370)

New England Biolabs

## Abstract

**Citation:** New England Biolabs Mixture for PCR Amplification (E7370). [protocols.io](https://doi.org/10.17504/protocols.io.cg7tzm)

[dx.doi.org/10.17504/protocols.io.cg7tzm](https://doi.org/10.17504/protocols.io.cg7tzm)

**Published:** 28 May 2015

## Protocol

### Step 1.

Adaptor Ligated DNA Fragments 23 µl


### Step 2.

NEBNext High Fidelity 2X PCR Master Mix 25 µl

 [AMOUNT](#)

25 µl Additional info:

 [REAGENTS](#)

 NEBNext High-Fidelity 2X PCR Master Mix - 50 rxns [M0541S](#) by [New England Biolabs](#)

### Step 3.

Index Primer 1 µl

### Step 4.

Universal PCR Primer 1 µl