

Loading Marker Mixture

New England Biolabs

Abstract

Citation: New England Biolabs Loading Marker Mixture. [protocols.io](https://www.protocols.io)

[dx.doi.org/10.17504/protocols.io.cq5vy5](https://doi.org/10.17504/protocols.io.cq5vy5)

Published: 03 Feb 2015

Protocol

Step 1.

Distilled water, **4 μ l**

 [AMOUNT](#)

4 μ l Additional info:


Step 2.

6X Purple no-SDS Loading Dye, **1 μ l**

 [AMOUNT](#)

1 μ l Additional info:

 [REAGENTS](#)

 Gel Loading Dye, Purple (6X), no SDS - 4.0 ml [B7025S](#) by [New England Biolabs](#)

Step 3.

DNA Ladder, **1 μ l**

 [AMOUNT](#)

1 μ l Additional info: