



Oct 12, 2018

Working

METADAP genetic protocol [↗](#)

PLOS One

[eric deflesselle](#)¹¹Université Paris Sud (Paris XI)[dx.doi.org/10.17504/protocols.io.ti2ekge](https://doi.org/10.17504/protocols.io.ti2ekge) [eric deflesselle](#)

ABSTRACT

Genomic DNA was extracted from circulating blood leucocytes by using Gentra Puregene Blood Kits according to the manufacturer's protocol (Qiagen) and was stored at -20°C.

The TRKB polymorphisms were genotyped by the IntegraGen company (Evry, France) using the Fluidigm® BioMark™ MHD system (Fluidigm Corporation, CA, USA), using TaqMan allelic discrimination (Coulbault et al., 2006), with the ABI Prism® 7900HT Sequence Detection System (Life Technologies).

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0205648>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Deflesselle E, Colle R, Rigal L, David DJ, Vievard A, Martin S, Becquemont L, Verstuyft C, Corruble E (2018) The TRKB rs2289656 genetic polymorphism is associated with acute suicide attempts in depressed patients: A transversal case control study. PLoS ONE 13(10): e0205648. doi: [10.1371/journal.pone.0205648](https://doi.org/10.1371/journal.pone.0205648)

PROTOCOL STATUS

Working

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited