

# Co-cultivation protocol for transforming wheat stripe rust

Yiheng Hu

## Abstract

Protocol modified for *Puccinia striiformis* f. sp. *tritici* (Pst), based on:

Michielse, C. B., Hooykaas, P. J., van den Hondel, C. A., & Ram, A. F. (2008). Agrobacterium-mediated transformation of the filamentous fungus *Aspergillus awamori*. *Nature protocols*, 3(10), 1671-1678.

I am grateful for critical suggestions from Prof. John Rathjen and Dr. Benjamin Schwessinger.

**Citation:** Yiheng Hu Co-cultivation protocol for transforming wheat stripe rust. **protocols.io**

dx.doi.org/10.17504/protocols.io.f9tbr6n

**Published:** 13 Mar 2017

## Guidelines

### Chemical solutions:

All chemicals used were from Sigma-Aldrich of the highest obtainable grades.

Buffers are best when fresh and not older than 3-6 months

- 0.2 M Acetosyringone (AS): Dissolve 785 mg of AS in DMSO to make up a total volume of 20 ml; aliquot and store in the dark at -20 °C. Do not thaw and use an aliquot more than twice, as activity of AS decreases during thawing and freezing
- 1% CaCl<sub>2</sub> (wt/vol): Dissolve 10 g of CaCl<sub>2</sub>•2H<sub>2</sub>O in water to make up a total volume of 1 liter; autoclave.
- 0.01% FeSO<sub>4</sub> (wt/vol): Dissolve 0.1 g of FeSO<sub>4</sub>•7H<sub>2</sub>O in water to make up a total volume of 1 liter; filter-sterilize.
- 20% glucose (wt/vol): Dissolve 200 g of glucose in water to make up a total volume of 1 liter; autoclave.
- 50% glycerol: Add 50 ml of glycerol to 50 ml of water to make up a total volume of 100 ml; autoclave.
- 1.25 M KH<sub>2</sub>PO<sub>4</sub>: Dissolve 170.1 g of KH<sub>2</sub>PO<sub>4</sub> in water to make up a total volume of 1 liter; autoclave.
- 1.25 M K<sub>2</sub>HPO<sub>4</sub>: Dissolve 217.7 g of K<sub>2</sub>HPO<sub>4</sub> in water to make up a total volume of 1 liter; autoclave.
- 1 M MES: Dissolve 195.24 g of MES in water to make up a total volume of 1 liter; adjust to pH

5.5 by adding NaOH; filter-sterilize. Solution can be stored for a month in the dark or alternatively aliquot and frozen at -20°C.

- K-buffer: Add solution 1.25 M  $\text{KH}_2\text{PO}_4$  to 1.25 M  $\text{K}_2\text{HPO}_4$  until pH 4.8 is reached.
- MN buffer: Dissolve 30 g of  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  and 15 g of NaCl in water to make up a total volume of 1 liter; autoclave.
- 20%  $\text{NH}_4\text{NO}_3$  (wt/vol): Dissolve 200 g of  $\text{NH}_4\text{NO}_3$  in water to make up a total volume of 1 liter; autoclave.
- Trace elements for IM medium: Dissolve 100 mg of  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ , 100 mg of  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ , 100 mg of  $\text{H}_3\text{BO}_3$ , 100 mg of  $\text{MnSO}_4 \cdot \text{H}_2\text{O}$  and 100 mg of  $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$  in water to make up a total volume of 1 liter; autoclave.

### Reagent set up:

- Hybond N+ filter (Amersham) need to be cut into small stripes at 5 mm in width. Do not autoclave.
- Induction medium (IM) (liquid): Add 0.8 ml of K-buffer, 20 ml of MN buffer, 1 ml of 1% (wt/vol)  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ , 10 ml of 0.01% (wt/vol)  $\text{FeSO}_4$ , 5 ml of trace elements for IM medium, 2.5 ml of 20% (wt/vol)  $\text{NH}_4\text{NO}_3$ , 10 ml of 50% (vol/vol) glycerol, 40 ml of 1 M MES, pH 5.5, and 10 ml of 20% (wt/vol) glucose to 900.7 ml of sterilized water to make up 1 liter of liquid IM. No glucose IM just substitute 10 ml 20% glucose by 10 ml  $\text{H}_2\text{O}$ .
- IM (solid): Dissolve 15 g of bacto agar in water to make up a total volume of 905.7 ml; autoclave. Add 0.8 ml of K-buffer, 20 ml of MN buffer, 1 ml of 1% (wt/vol)  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ , 10 ml of 0.01% (wt/vol)  $\text{FeSO}_4$ , 5 ml of trace elements for IM medium, 2.5 ml of 20% (wt/vol)  $\text{NH}_4\text{NO}_3$ , 10 ml of 50% (vol/vol) glycerol, 40 ml of 1 M MES, pH 5.5, and 5 ml of 20% (wt/vol) glucose to make up 1 liter of solid IM.

Wheat cultivar used: *Morocco*

*A. tumefaciens* strain: GV3101 pM90

*Pst* strain used: *P. striiformis* f.sp. *tritici* strain 104E137A-

### *Pst* spore reproduction and collection:

Spores are resuspended in 3M™ Novec™ 7100 engineered fluids immediately before inoculation. 7-day old seedlings were inoculated with fresh urediniospores of *Pst* using flat brush and incubated for 24 h in 100% humidity at 8°C in the dark. Subsequently plants were transferred to a growth chamber at 17°C with a 16:8 light cycle. Spores can be harvested after 14 dai, with 2 days intervals, by dusting the heavily infected leaves on aluminium foil or a weighing paper.

### GUS staining of *Pst* spores:

substrate solution: 0.05M sodium phosphate buffer, pH7.4; 1mM EDTA; 0.1% (v/v) TritonX-100; 0.1mM K<sub>3</sub>Fe(CN)<sub>6</sub>; 0.5mg X-Gluc/ml

Freshly harvested *Pst* spores are germinated on the surface of water for 6 h, before replacing all water by the substrate solution using glass Pasteur pipette. Spores in the substrate solution are then incubated at 37 °C for 40 h, before moving to a glass slide using glass Pasteur pipette for imaging.

**DNA extraction from *Pst* spores** is according to Yiheng Hu: High quality DNA extraction from Fungi\_small scale. protocols.io  
<https://www.protocols.io/view/high-quality-dna-extraction-from-fungi-small-scale-exmbfk6>

### Critical steps to obtain high transformation efficiency:

- Must use fresh spores/wheat seedlings/*Acetosyringone/Agrobacterium*
- Must use glass plates/beakers to germinate spores
- Make sure the filter contain co-cultured spores attached tightly to the wheat leaf to maximize the infection.

### Before start

1. Grow 12 pots of wheat seedlings (15 seeds/pot) and infect them with 50 µL of *Pst* spores to reproduce fresh spores for co-culture.
2. Plant certain numbers of wheat (depend on each experiment, I do 6 to 8 pots) 7 days before the experiment start (day 1)
3. Streak out *A. tumefaciens* contains the T-DNA in binary vector onto a LB-Agar plate and grow at 28°C for 48 h. Pick up one colony of the *A. tumefaciens* into 5 ml liquid LB media with proper selection and grow overnight. The co-culture experiment need to start from a nice overnight *A. tumefaciens* culture.

### Protocol

#### Step 1.

Centrifuge the overnight *A. tumefaciens* culture at 5000xg for 10 mins, wash the pellet using induction media (IM, with glucose) once and resuspend the pellet with 5 mL IM (with glucose).

 **DURATION**

00:10:00

#### Step 2.

Add 2ml of Agro-IM solution into 8ml liquid IM (with glucose) in a 50 mL falcon tube and put it at 28°C, shake until OD<sub>600</sub>=0.6 (6 h depends on *Agrobacterium* strain)

### Step 3.

Harvest fresh spores from previous infected wheat, estimate the number of spores based on a counting experiment:

I first suspended them 0.1% (v/v) of the non-ionic detergent Triton X-100 to wet and separate them. Then I diluted the spore solution 10, 100 and 1000 times, and counted *Pst* spores using a haemocytometer I also estimated the spore volume by comparing the volume of a known number of spores with equivalent volumes of water in an Eppendorf tube. I estimated that 10 mg (20 ml) of *Pst* spores equals  $3.5 \times 10^6$  spores.

### Step 4.

Dust 100ul of spores onto a glass plate/beaker at 14 cm in diameter (indicated by the image below), gently spray water on to the surface using a water sprayer to help the spores distribute equally. Cover the plate using parafilm and put it into 8°C for 6 h




### Step 5.

After the IM solution containing *A. tumefaciens* reaches OD=0.6, centrifuge at room temperature at 5000xg for 10mins, resuspend it using 10 ml liquid IM (without glucose) with 200 uM Acetosyringone (AS) added.


### Step 6.

Take out four IM-Agar plates (with 200uM AS added) and take out the spores from 8°C. Using sterilized tweezers to tweeze the cut hybond N+ filter to attach spores from the surface. Dry the filter a little using the airflow in the hood for 1min, and put it onto IM-Agar (spore side up). Normally five pieces of filter strips per plate and two plates per pot of wheat seedlings ( 15 seedlings/pot).



### Step 7.

Pipette 100ul IM (no glucose, AS added) containing *A. tumefaciens* onto each filter with spores, wait them to dry for 3mins and put them at 20°C to process the co-culture.




 **DURATION**

36:00:00

### Step 8.

After co-cultivation of 36 h, take out the plates from the constant temperature room and prepare previously grown 6-8 pots of fresh 9-day old wheat seedlings. Tweeze the filter up from plate, gently spray water onto it and attach onto the wheat leaf (primary leaf). After all filter is attached, gently spray a little water onto the plants and then put them into 100% humidity (a container with water) at 8°C to germinate 24 h again.



#### DURATION

24:00:00

#### **Step 9.**

Transfer the infected wheat into growth chamber (17°C 16:8 light:dark cycle, 60% humidity).

#### **Step 10.**

After 8 days, check the infection of spores on the wheat, and take off the filter papers to let pustules erupt.

Plant 8 pots of wheat to amplify the transformed spores for selection

#### **Step 11.**

Collect spores, reinfect a small number of spores for amplification and further analysis (GUS staining, DNA extraction etc.).