

CviRI Buffer B

DAVID DUNIGAN AND IRINA AGARKOVA

Abstract

For use in [CviRI Purification From XZ-6E Virus Infected NC64A Chlorella](#).

Citation: DAVID DUNIGAN AND IRINA AGARKOVA CviRI Buffer B. **protocols.io**

dx.doi.org/10.17504/protocols.io.esxbefn

Published: 13 Jun 2016

Protocol

Step 1.

20 mM Tris-HCl, pH 7.8

Step 2.

0.5 mM EDTA

Step 3.

7 mM 2-ME, 10% Glycerol