

# 10X MOPS Buffer

DAVID DUNIGAN AND IRINA AGARKOVA

## Abstract

For use in [Preparation of Denaturing Agarose Gel for RNA Analysis](#).

**Citation:** DAVID DUNIGAN AND IRINA AGARKOVA 10X MOPS Buffer. **protocols.io**

[dx.doi.org/10.17504/protocols.io.etgbejw](https://dx.doi.org/10.17504/protocols.io.etgbejw)

**Published:** 13 Jun 2016

## Protocol

### Step 1.

10X MOPS consists of the following (per liter):

### Step 2.

0.2 M MOPS pH 7 with NaOH

### Step 3.

50 mM sodium acetate

### Step 4.

10 mM EDTA