

SYBR Gold Working Stock

B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John

Abstract

For use in [Wet-mount Method for Enumeration of Aquatic Viruses](#)

Citation: B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John SYBR Gold Working Stock. **protocols.io**
dx.doi.org/10.17504/protocols.io.c82zyd

Published: 07 Jan 2016

Protocol

Step 1.

Dilute SYBR Gold (Invitrogen; 10,000X stock) into PBS (phosphate buffered saline) to prepare 1000x solution. Prepare fresh daily. Store in dark between uses.