



Oct 06, 2019

## PCR and cloning

Norfitriah Mohamed Sohaimi<sup>1</sup>, Mohd Hair Bejo<sup>1</sup>, Abdul Rahman Omar<sup>1</sup>, Aini Ideris<sup>1</sup>, Nurulfiza Mat Isa<sup>2</sup>

<sup>1</sup>Faculty of Veterinary Medicine, Universiti Putra Malaysia, Serdang, Selangor, Malaysia., <sup>2</sup>Faculty of Biotechnology and Biomolecular Sciences, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.

1

Works for me

dx.doi.org/10.17504/protocols.io.7zfh3n



Norfitriah Mohamed Sohaimi



### MATERIALS

NAME

CATALOG #

VENDOR

TOPO<sup>®</sup>; TA Cloning<sup>®</sup>; Kit, with pCR<sup>®</sup>; 2.1-TOPO<sup>®</sup>; One Shot TOP10 Chemically Competent *E. coli*, and PureLink<sup>®</sup>; Quick Plasmid Miniprep Kit

K450002

Thermo Fisher

### MATERIALS TEXT

Amplification of DNA was performed by two set of published hexon gene primers, H1/H2 and H3/H4 [Raue and Hess, 1998] and fiber gene primer, FibF/FibR, FibF: 5'-GGTCTACCCCTTTGGCTCC-3' and FibR: 5'-GCGTCGTAGATGAAGGGAGG-3' [Norfitriah et al., 2018] according to manufacture protocol (Bioline, UK). The PCR products were analyzed by electrophoresis in a 1% agarose gel stained with RedSafe<sup>™</sup> Nuclei Acid Staining solution (iNtRON, Korea) at 70 volts for 45 minutes and visualized under U.V. transillumination. Purification of PCR products were carried out by using MEGAquick-spin<sup>™</sup> Total Fragment DNA Purification kit (iNtRON) based on the manufacture recommendation. Purified PCR products from H1/H2 and FibF/FibR were cloned into the pCR<sup>™</sup> 2.1-TOPO<sup>®</sup> vector using TOPO TA Cloning kit (Invitrogen, USA). Positive cloned was analyzed by colony-PCR prior plasmid extraction using DNA-spin<sup>™</sup> Plasmid DNA Purification kit (iNtRON) and stored at -20°C until used.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited