



Phage Buffer

Dr. Steven Wilhelm¹

¹Microbiology Department. University of Tennessee, Knoxville.

[dx.doi.org/10.17504/protocols.io.iitccen](https://doi.org/10.17504/protocols.io.iitccen)

The Aquatic Microbial Ecology Research Group - AMERG (The Buchan, Zinser and Wilhelm labs)

CyanoHABs



Ashley Humphrey

University of Tennessee, Knoxv...



ABSTRACT

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

SAFETY WARNINGS

See SDS (Safety Data Sheet) for safety warnings and hazards.

- 1 In a one liter bottle, add 10 ml MgSO_4 .

 10 ml 1 Molar

- 2 Add 1 ml CaCl_2

 1 ml 1 Molar

- 3 Add 1 ml Tris-HCl (pH~7)

 1 ml 1 Molar

- 4 Add 1 g gelatin.

 1 g

- 5 Add Milli-Q water to a final volume of 1000 mL and autoclave sterilize.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited