

Stellaris® RNA FISH Fixation Buffer

LGC Biosearch Technologies

Abstract

Citation: LGC Biosearch Technologies Stellaris® RNA FISH Fixation Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.cunwvd

Published: 23 Feb 2016

Protocol

Step 1.

37% Formaldehyde solution

 **AMOUNT**

1 ml Additional info:

Step 2.

10X Phosphate Buffered Saline (PBS), RNase-free

 **AMOUNT**

1 ml Additional info:

Step 3.

Nuclease-free water

 **AMOUNT**

8 ml Additional info: