

Euplotes crassus transformation through microinjection into the macronucleus

Rachele Cesaroni

Abstract

Citation: Rachele Cesaroni Euplotes crassus transformation through microinjection into the macronucleus. protocols.io

dx.doi.org/10.17504/protocols.io.g3fbyjn

Published: 30 Jan 2017

Protocol

Step 1.

Euplotes crassus cells of two different mating types were diluted 1:10 in artificial sea water and were fed with E.coli (1:1000).

📌 NOTES

Rachele Cesaroni 26 Jan 2017

Recipe for complete seawater (1 L):

36 g Reef Crystals

1 ml Walne's solution

1 ml of 10 µg/ml FeSO₄

0.2 ml of 2 mg/ml thiamine (light sensitive; store at 4°C)

Add distilled water to 1 L

Step 2.

Cells of both mating types were grown at 24°C for 4 days with a 12h light/12h dark cycle and then mixed at 24°C (after few hours they started to form pairs).

■ ANNOTATIONS

Rachele Cesaroni 30 Jan 2017

Aeration was provided while cells were growing, but not while cells were mating.

Step 3.

50h after formation of the pairs (donut cells) cells were isolated and put into 2% BSA diluted in artificial sea water in order to prepare drops for microinjection (ideally one cell each drop).

📌 NOTES

Rachele Cesaroni 30 Jan 2017

Cells have a donut shape (anlagen) and macronucleus is easier to inject.
50h after formation of the pairs a round of amplification of the genome occurs.

Step 4.

When drops were ready a layer of Mineral Oil was put above to not let the drops evaporate.

📌 NOTES

Rachele Cesaroni 30 Jan 2017

Cells have a donut shape (anlagen) and macronucleus is easier to inject.
50h after formation of the pairs a round of amplification of the genome occurs.

Step 5.

Phenol/chloroform extracted, precipitated and resuspended in injection buffer (114 mM KCl, 20 mM NaCl, 3 mM sodium phosphate buffer) DNA (artificial nanochromosome with telomeres, but without 3' overhangs) 3-5 µg/µl was injected into the macronucleus using Eppendorf Femtotips II injection needle.

Step 6.

Cells were then recovered in 400 µl of sea water plus E.coli (0.4 µl) at put at 24°C.