

hMSC basal media

Abstract

hMSC Basal Media

Citation: hMSC basal media. [protocols.io](https://www.protocols.io)

[dx.doi.org/10.17504/protocols.io.ccrsv5](https://doi.org/10.17504/protocols.io.ccrsv5)

Published: 18 Jun 2014

Protocol

Step 1.

86% DMEM

 [AMOUNT](#)

172 ml Additional info:

Step 2.

1% Pen Strep

 [AMOUNT](#)

1 ml Additional info:

Step 3.

1% NEAA

 [AMOUNT](#)

2 ml Additional info:

Step 4.

2% L-Glutamine

 [AMOUNT](#)

4 ml Additional info:

Step 5.

10% FBS

 [AMOUNT](#)

20 ml Additional info: