CITE-seq Protocols

Marlon Stoeckius, Peter Smibert

Abstract


This collection contains our main protocols for performing CITE-seq and Cell Hashing, specifically on <u>Drop-seq</u> or <u>10x Genomics single cell 3' v2 chemistry</u>.

CITE-seq:

Cellular Indexing of Transcriptomes and Epitopes by Sequencing (<u>CITE-seq</u>) is a multimodal single cell phenotyping method developed in the <u>Technology Innovation lab</u> at the New York Genome Center in collaboration with the Satija lab.



CITE-seq uses DNA-barcoded antibodies to convert detection of proteins into a quantitative, sequenceable readout. Antibody-bound oligos act as synthetic transcripts that are captured during most large-scale oligodT-based scRNA-seq library preparation protocols (e.g. 10x Genomics, Drop-seq, ddSeq).




This allows for immunophenotyping of cells with a potentially limitless number of markers and unbiased transcriptome analysis using existing single-cell sequencing approaches.

Cell Hashing:

Sample multiplexing and super-loading on single cell RNA-sequencing platforms.

<u>Cell Hashing</u> uses a series of oligo-tagged antibodies against ubiquitously expressed surface proteins with different barcodes to uniquely label cells from distinct samples, which can be subsequently pooled in one scRNA-seq run. By sequencing these tags alongside the cellular transcriptome, we can assign each cell to its sample of origin, and robustly identify doublets originating from multiple samples.



Citation: Marlon Stoeckius, Peter Smibert CITE-seq Protocols. protocols.io

dx.doi.org/10.17504/protocols.io.ngzdbx6

Published: 29 May 2018

Collection

PROTOCOLS

1. CITE-seq

CONTACT: Peter Smibert

2. Cell Hashing

CONTACT: Peter Smibert

3. CITE-seq and Cell Hashing

CONTACT: Peter Smibert