

PBS / DAPI 1:2000 Dilution Mixture

Alan Cone

Abstract

Citation: Alan Cone PBS / DAPI 1:2000 Dilution Mixture. **protocols.io**

dx.doi.org/10.17504/protocols.io.eifbcbn

Published: 10 Feb 2016

Protocol

Step 1.

Add 1 mL 1 x PBS to a 1.5 mL microcentrifuge tube.

AMOUNT

1 ml Additional info:

REAGENTS

✓ 1X PBS (Phosphate-buffered saline) by Contributed by users

Step 2.

Add 0.5 µL of a 2.5 mg/mL (or a 1:2000 dilution) of DAPI to the 1 x PBS.

REAGENTS

✓ DAPI (2.5mg/mL) by Contributed by users