

Introduction to read mapping

Frank Aylward

Abstract

This is an example of a simple read mapping workflow. It is designed to be performed via the command line on an Ubuntu 16.06 OS.

After completing this tutorial you should:

- 1) Have a practical understanding of how read mapping analyses are performed in the command line.
- 2) Understand the basics of how to process SAM and BAM files.
- 3) Be able to calculate depth of coverage of a contig/scaffold/chromosome in a query read dataset.

Citation: Frank Aylward Introduction to read mapping. **protocols.io**

dx.doi.org/10.17504/protocols.io.pgwdjxe

Published: 17 Apr 2018

Protocol

Ensure the appropriate tools are installed

Step 1.

First ensure that the following tools are installed and are in your PATH:

Bowtie 2 version 2.2.6

samtools Version: 0.1.19-96b5f2294a

bedtools v2.25.0

The exact versions may not be critical, but those used in designing this tutorial are listed for completeness.

First we need to get a reference genome. In this case we will be working with a bacterium called

Marinimicrobial UBA

Step 2.

First we need to get a reference genome to map reads against. We'll download the genome of Marinimicrobia UBA2153 here.

Download using the unix command wget:

```
wget  
ftp://ftp.ncbi.nlm.nih.gov/genomes/all/GCA/002/328/885/GCA_002328885.1_ASM232888v1/GCA_002328885.1_ASM232888v1_genomic.fna.gz
```

And since the fasta file is compressed we will use the Unix tool gunzip to decompress it:
gunzip GCA_002328885.1_ASM232888v1_genomic.fna.gz

Build bowtie2 reference

Step 3.

To use the read mapping tool bowtie2 we will need to first index the fasta file. This can be done using the command bowtie2-build. It will create several index files with different suffixes and the prefix that we give in the command after the fasta file (in this case UBA2153).

```
bowtie2-build GCA_002328885.1_ASM232888v1_genomic.fna UBA2153
```

Get the reads for mapping

Step 4.

Now we need to get the reads that we will use for mapping.

```
fastq-dump -X 1000000 --gzip --split-3 SRR5322088
```

Map the reads with bowtie2

Step 5.

```
bowtie2 -1 SRR5322088_1.fastq -2 SRR5322088_2.fastq -x UBA2153 -S mapping_output.SAM
```

Now process the SAM file created by bowtie2 with samtools

Step 6.

```
samtools view -bS -F 4 mapping_output.SAM > mapping_output.bam  
samtools sort mapping_output.bam mapping_output.sort  
samtools index mapping_output.sort.bam
```

samtools idxstats mapping_output.sort.bam