



May 10,
2019

Working

U Mass - Body composition (organs) [↗](#)

Jason Kim¹

¹University of Massachusetts

dx.doi.org/10.17504/protocols.io.xsnfnde

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

The EchoMRI 3-in-1 uses ¹H- magnetic resonance spectroscopy to noninvasively measure fat mass, lean mass, and water mass in individual organs. Fat, lean, and water mass composition in specific organs is altered in obesity.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=148&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
EchoMRI 3-in-1		Echo Medical Systems		
Pentobarbital	NDC76478-501-50	Oak Pharmaceuticals, Inc.		

MATERIALS TEXT

Note:

[EchoMRI](#), [RRID:SCR_017104](#)

- 1 Euthanize mice using intraperitoneal injection of pentobarbital.
- 2 Quickly dissect and extract specific organs to measure their composition.
- 3 Place specific organ into an instrument column.
- 4 Measure composition of organ using instrument standard operating procedure.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited