

# Ascorbic Acid Antifade Solution

**B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John**

## Abstract

For use in [Wet-mount Method for Enumeration of Aquatic Viruses](#).

**Citation:** B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John Ascorbic Acid Antifade Solution.  
**protocols.io**

[dx.doi.org/10.17504/protocols.io.c87zzm](https://dx.doi.org/10.17504/protocols.io.c87zzm)

**Published:** 07 Jan 2016

## Guidelines

This solution is needed only if analyzing unconcentrated samples.

## Protocol

### Step 1.

Dissolve 0.1 g ascorbic acid into 1 mL PBS creating a 10% (wt/vol) solution. Prepare fresh daily.