

# TKN Buffer

Christopher M. Bartley

## Abstract

This buffer was modified from Lindsay and McCaffery 2014 (DOI: 10.1111/boc.201200076). We use to investigate FMRP protein-protein interactions.

**Citation:** Christopher M. Bartley TKN Buffer. **protocols.io**

[dx.doi.org/10.17504/protocols.io.cvhw35](https://dx.doi.org/10.17504/protocols.io.cvhw35)

**Published:** 28 Mar 2015

## Protocol

### Step 1.

20mM Tris pH 7.5

### Step 2.

100mM KCl

### Step 3.

0.5% Triton X-100

### Step 4.

100uM Okadaic Acid to if investigating phosphorylation-dependent FMRP protein-protein interactions

### Step 5.

1X HALT

### Step 6.

8U/10mL DNase 1