T4 PNK Minus Reaction

New England Biolabs

Abstract

Citation: New England Biolabs T4 PNK Minus Reaction. protocols.io

dx.doi.org/10.17504/protocols.io.cp3vqm

Published: 31 Jan 2015

Protocol

Step 1.

DNA up to 300pmol of 5' termini

Step 2.

10X T4PNK Reaction Buffer, **5 μl**

■ AMOUNT

5 μl Additional info:

Step 3.

10mM ATP, **5 μl**

■ AMOUNT

5 μl Additional info:

Step 4.

T4 PNK (3' phosphatase minus), **1 μl** (10 units)

■ AMOUNT

1 μl Additional info:

Step 5.

Water up to 50 µl