Published: 31 Jan 2015

Radioactive T4PNK Minus Reaction

New England Biolabs

Abstract

Citation: New England Biolabs Radioactive T4PNK Minus Reaction. protocols.io

dx.doi.org/10.17504/protocols.io.cqavsd

Published: 31 Jan 2015

Protocol

Step 1.

DNA, up to **50 pmol** of 5' termini

Step 2.

10X T4PNK Reaction Buffer, **5 μl**

■ AMOUNT

5 μl Additional info:

Step 3.

ATP **50 pmol** of [γ -32P] ATP, **5** μ l

■ AMOUNT

5 μl Additional info:

Step 4.

T4PNK 3' phosphatase minus, 2 μl (20 units)

Step 5.

Water up to 50 µl