

1X CviJI Assay Buffer

DAVID DUNIGAN AND IRINA AGARKOVA

Abstract

For use in [CviJI Purification From IL-3A Virus Infected NC64A Chlorella](#).

Citation: DAVID DUNIGAN AND IRINA AGARKOVA 1X CviJI Assay Buffer. **protocols.io**
dx.doi.org/10.17504/protocols.io.espbcdn

Published: 13 Jun 2016

Protocol

Step 1.

20 mM GlycylGlycine, pH 8.5 (with KOH)

Step 2.

10 mM MgAc

Step 3.

0.1 mM DTT

Step 4.

50 mM KOAc

Step 5.

100 µg/mL casein (optional)